

Four-Step Analysis of Smokeless Tobacco

Tobacco has been studied extensively using analytical pyrolysis, frequently from the standpoint of identifying specific compounds in the smoke. Whether or not the tobacco product is intended to be smoked, however, significant information may be obtained using a series of thermal sampling steps. The tobacco shown here is a flavored, smokeless product, flavored with wintergreen.

At 100°C, two of the principal components may already be seen, as shown in Figure 1. Nicotine, naturally present, and methyl salicylate, wintergreen flavoring, an additive, are readily apparent. In addition, vitamin E, a natural product of tobacco, is also seen.

At 200°C, acetic acid and propylene glycol are detected. At this temperature, the tobacco begins to degrade, and pyrolysis products may be seen. At 300°, degradation products of cellulose become significant, including furans and levoglucosan.

There are several advantages in performing multiple thermal steps in this analysis. Some of the minor constituents (like the vitamin E) would be difficult to identify if a single pyrolysis run, producing all the products at once, were done instead. This approach also makes it clear that the nicotine is actually quite volatile and, even when tobacco is smoke, is not a pyrolysis product, but simply volatilized. Additives, such as menthol, methyl salicylate, propylene glycol and glycerine, may also be determined easily, apart from the products generated when Tobacco is pyrolyzed.




Figure 1. Analysis of wintergreen flavored smokeless tobacco at 100, 200, 300 and 600°C.

Table I Peak Identification

1. Methyl salicylate
2. Nicotine
3. Vitamin E
4. Acetic acid
5. Propylene glycol
6. 2-Furanmethanol
7. Hydroxymethyl furancarboxaldehyde
8. Levoglucosan
9. Toluene
10. Phenol

Equipment

The tobacco sample was analyzed directly from the container, using a CDS Model 5200 Pyroprobe, interfaced to a Clarus 500 gas chromatograph/mass spectrometer.

Model 5200 Conditions

Valve Oven: 300°C
Transfer Line: 325°C
Temperature: As indicated
Time: 15 seconds

The 5200 was used in the direct pyrolysis mode, without trapping.

GC Conditions

Carrier: Helium
Column: Rxi-5ms (30m X 0.25mm)
Detector: Clarus 500 MS

GC Program:

Initial: 40°C for 2 minutes
Ramp: 10°C/min.
Final: 300°C

FOR MORE INFORMATION
CONCERNING THIS APPLICATION,
WE RECOMMEND THE
FOLLOWING READING:

E. B. Sanders, et al., A model that distinguishes the pyrolysis of D-glucose, D-fructose, and sucrose from that of cellulose. Application to the understanding of cigarette smoke formation, *J. Anal. Appl. Pyrolysis*, 66 (2002) 29-50.

Additional literature on this and related applications may be obtained by contacting your local CDS Analytical representative, or directly from CDS at the address below.

CDS Analytical, Inc. has been a leader in the design and manufacture of laboratory instruments for sample preparation and analysis since 1969. We are dedicated to providing the best possible instruments for both research and routine analysis. Well known in the field of pyrolysis, CDS manufactures the Pyroprobe® 5000, 5150, 5200 and 5250 autosampler for the introduction and analysis of solid materials by GC, MS and FT-IR. CDS offers a complete line of dynamic headspace instruments for the analysis of volatile organic compounds in environmental, pharmaceutical and food applications, including the model 8400 four-position autosampler for complex, multicomponent materials investigation. Our customers, their requirements and applications are important to us. To help meet your needs, we offer a wide range of analytical information and the services of our applications laboratory. If you would like additional information, please contact us at the address below, call us at 1 800 541 6593, or log onto cdsanalytical.com.