

CDSolutions


APPLICATIONS INFORMATION USING ADVANCED SAMPLE HANDLING TECHNOLOGY

The Use of Pyrolysis GC/MS to Characterize Indoor Dust

As indoor dust may contain toxic materials along with allergen causing natural polymers from spiders, insects and mites, the study of indoor dust can be used to determine potential environmental health risks in the environment. Using a CDS Pyroprobe as an injection strategy for gas chromatography/mass spectrometry, we are able to understand both the semi-volatile and polymer nature of dust.

We took dust from kitchens, bathrooms, and bedrooms in homes, and ran them on the GC/MS using a CDS Pyroprobe 5200 to introduce the samples. Small amounts of dust (<100µg) were desorbed at 300°C, and then pyrolyzed at 750°C. As expected, dust collected turned out to be a complex mixture. In this study, just a few components were identified. Each location had pyrolysis products such as phenols from wool and hair, levoglucosan and furans from cellulose (cotton or paper), benzoic acid from polyester (clothing and carpeting), and styrene from polystyrene. Furthermore, not all locations were identical; each contained different types and relative amounts of products, relating to the items found in the area. For example, home "K", furnished with polyester carpeting, had more polyester; and home "T" had more cellulose from cotton clothing.

Semivolatiles present included fatty acids and cholesterol from skin and cooking oils, and the plasticizer bis(2-ethylhexyl) phthalate. By sampling dust, and items suspected of contributing to dust, thermal desorption and pyrolysis GC/MS can be used to help identify dust make-up.



Peaks from Wool or Hair(W), Cotton (C), and Polyester (P) labeled in four dust samples.

Equipment

CDS Model 5200 Pyroprobe in direct (no trapping) mode interfaced to a Varian 3800 GC Saturn 2000 Ion Trap detector.

Model 5200 Conditions

Desorption:

Valve Oven: 325°C

Transfer Line: 325°C

Pyroprobe Temperature: 300°C

Pyroprobe Time: 15 seconds

Interface Final: 300°C for 3 minutes

Pyrolysis:

Valve Oven: 325°C

Transfer Line: 325°C

Pyroprobe Temperature: 750°C

Pyroprobe Time: 15 seconds

Interface Final: 325°C for 3 minutes

GC Conditions

Carrier: Helium

Injector: 325°C

Split: 50:1

Column: 5% phenyl (30m X 0.25mm)

Detector: Ion Trap

Range: 35 - 550amu

GC Program:

Initial: 40°C for 2 minutes

Ramp: 10°C/min.

Final: 300°C for 5 minutes

FOR MORE INFORMATION
CONCERNING THIS APPLICATION,
WE RECOMMEND THE
FOLLOWING READING:

K. Jansson et. al., Multistep thermal characterization of polymers using GC-MS, American Lab, 39 (6) (2007) 16-19.

Additional literature on this and related applications may be obtained by contacting your local CDS Analytical representative, or directly from CDS at the address below.

CDS Analytical, Inc. has been a leader in the design and manufacture of laboratory instruments for sample preparation and analysis since 1969. We are dedicated to providing the best possible instruments for both research and routine analysis. Well known in the field of pyrolysis, CDS manufactures the Pyroprobe® 5000, 5150, 5200 and 5250 autosampler for the introduction and analysis of solid materials by GC, MS and FT-IR. CDS offers a complete line of dynamic headspace instruments for the analysis of volatile organic compounds in environmental, pharmaceutical and food applications, including the model 8400 four-position autosampler. CDS also manufactures the Dynatherm line of thermal desorption instruments including the 9000 series for air monitoring and the 9300 TDA. Our customers, their requirements and applications are important to us. To help meet your needs, we offer a wide range of analytical information and the services of our applications laboratory. If you would like additional information, please contact us at the address below, call us at 1 800 541 6593, or log onto **www. cdsanalytical.com**.