

Sensitive, High-Throughput *In Vitro* ADME P-Glycoprotein Inhibition With Agilent RapidFire/MS Systems

Application Note

Authors

Vaughn Miller
Agilent Technologies, Inc.
Wakefield, MA, USA

Introduction

Interactions between drug candidates and transporters are routinely evaluated during drug discovery. P-glycoprotein (P-gp) is a transporter which is commonly assessed at an early phase of the process to determine if the drug of interest is a substrate or an inhibitor of P-gp, so that poor drug candidates can be eliminated. Since analysis early in the drug discovery process involves large sample sets, a high-throughput approach is desirable. Agilent RapidFire/MS systems enable high-throughput sample processing to streamline ADME assay analysis.



Agilent Technologies

Using RapidFire High-Throughput Mass Spectrometry to Analyze P-Glycoprotein Inhibition

The RapidFire High-throughput Mass Spectrometry System was employed to assess the effect of the drug candidates on bidirectional transport of the substrate, digoxin, in a Caco-2 cell-based system which expresses the P-gp transporter (Figure 1).

RapidFire analysis of this assay achieved sample cycle times of six to ten seconds per injection, enabling an ultra-fast, bioanalytical method to assess P-gp inhibition.

Solid phase extraction (SPE) based RapidFire does not require the sample preparation of chromatography-based systems, which makes the system a straightforward high-throughput solution. Combined with mass spectrometry, the system delivers the sensitivity required to analyze P-gp

inhibition assays and demonstrates excellent correlation to LC/MS/MS and radioactive methods (Figure 2). Side-by-side testing of cyclosporin A P-gp inhibition IC_{50} in all three assay methods yielded similar results as shown in Figure 3. Furthermore, the RapidFire/MS/MS total workflow was more than 10 times faster than the LC/MS/MS method (Table 1), but demonstrated similar sensitivity, selectivity, reproducibility, linearity, and robustness.




Figure 1. Schematic of the P-gp inhibition assay using a Caco-2 cell model. The basolateral to apical (B to A) direction of the bi-directional assay is shown here.

Table 1. Throughput comparison between the RapidFire/MS/MS and LC/MS/MS analyses of P-gp inhibition samples.

Time (h)	RapidFire/MS/MS	LC/MS/MS
Sample Preparation	0.75	0.75
Sample Analysis	1.50	24.0
Data Review	0.25	0.50
Total	2.50	25.3




Figure 2. Correlation plots comparing % P-gp inhibition of test compounds in RapidFire/MS/MS, LC/MS/MS, and radiolabel assays.




Figure 3. Side-by-side testing of cyclosporin A P-gp inhibition IC_{50} in all three assay methods.

Conclusions

The RapidFire High-throughput Mass Spectrometry System demonstrated a number of key benefits for the assessment of early drug candidates in P-Glycoprotein inhibition assays: rapid sample processing speeds, increased throughput and laboratory efficiency, and equivalent inhibition results as compared to conventional LC/MS methods. As a result, incorporation of RapidFire/MS systems into the *in vitro* ADME phase of the drug discovery process enables efficiency and productivity advances unrivaled by other technologies.

References

1. Wagner, A. *et al.* Ultrafast mass spectrometry based bioanalytical method for digoxin supporting an *in vitro* P-glycoprotein(P-gp) inhibition screen. *Rapid Commun Mass Spectrom*, **2011**, 25:1231-1240.

www.agilent.com/lifesciences/rapidfire

For research purposes only and not for use in diagnostic procedures. The information described here is intended for reference and research purposes only. Agilent Technologies offers no guarantee as to the quality or suitability of this data for your specific application.

Information, descriptions, and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc. 2011
Published in the USA, September 15, 2011
5990-9082EN



Agilent Technologies