



Speciation of Inorganic Arsenic in Baby Rice Cereals Using HPLC-ICP-MS

Application Note

Food Safety

Authors

Rima Juskelis and Jack Cappozzo
Institute for Food Safety and Health
Illinois Institute of Technology, Bedford
Park, IL USA

Jenny Nelson

Agilent Technologies, Inc.
5301 Stevens Creek Blvd. Santa Clara,
CA, USA

Abstract

A survey of baby rice cereals using ICP-MS and the method specified in FDA EAM: Section 4.11 detected all arsenic species in all of the cereals tested. Limit of Detection (LOD) values for As species ranged from approximately 0.9 ng/g to 1.8 ng/g, and Limit of Quantitation (LOQ) from 7 µg/kg to 14 µg/kg, and %RSDs ranged from 1.7 to 6.7%.



Agilent Technologies

Introduction

Rice is one of the main sources of inorganic arsenic (iAs) consumed in the world population's diet [1]. Arsenic occurs naturally in the environment and can also be present as the result of human activity, such as introduction to the soil through the use of arsenic-based pesticides prior to the 1970s. Rice plants are especially efficient at accumulating arsenic from their environment, because the flooded areas in which they are grown make it easier to take up arsenic compounds [2,3].

Inorganic arsenic is a known carcinogen, and chronic exposure to low levels of arsenic has been linked to increased risk of bladder, lung, and skin cancer, as well as Type 2 diabetes and cardiovascular disease [4]. Brown rice tends to contain more arsenic than white rice [1,5,6]. Rice samples from the United States have shown higher inorganic arsenic levels compared to other samples from around the world, and infant rice cereal sold in the US is generally made from American rice [2].

This application note describes use of the high performance liquid chromatography/inductively coupled plasma mass spectrometry (HPLC-ICP-MS) method specified in the FDA Elemental Analysis Manual (EAM): Section 4.11 for the specification of arsenic in baby rice cereals. The arsenic species are separated using isocratic anion-exchange HPLC, and the ICP-MS is used as an arsenic-specific detector monitoring m/z 75 for arsenic-containing chromatographic peaks. It is operated in helium collision cell mode to eliminate interference from possible coeluting chloride species. Limits of detection and quantitation (LOD and LOQ) were below 15 $\mu\text{g/g}$ for all arsenic species, with relative standard deviations (%RSDs) all below 7%. The method has been used to provide information on arsenic levels in rice and rice-based infant cereals necessary to help understand the health risk to infants [7].

Experimental

Chemicals and Standards

All solutions were prepared using deionized water with resistance > 18 $\text{M}\Omega\cdot\text{cm}$ with a Milli-Q system (Millipore).

| Chemical/Standard | Source |
|--|----------------------|
| Nitric acid, OPTIMA ultra-pure grade | Fisher Scientific |
| Hydrogen peroxide, OPTIMA ultra-pure grade | Fisher Scientific |
| Ammonium phosphate | ACROS |
| Ammonium hydroxide, ultrapure | Fisher Scientific |
| Isopropanol, HPLC-grade | Fisher Scientific |
| Arsenite Stock Standard (AsIII), 1,000 mg/L | Spex Certiprep |
| Arsenate Stock Standard (AsV), 1,000 mg/L | Spex Certiprep |
| Monosodium acid methane arsonate (MMA, 98.5% purity) | Chem Service |
| Disodium acid methane arsonate (DMA, 98.9% purity) | Chem Service |
| Arsenobetaine (AsB) | Fluka |
| Rhodium/germanium internal standard | Agilent Technologies |
| ICP-MS Arsenic Standard | Plasma |

Table 1. Rice Flour Standard Reference Materials

| SRM/CRM | Total As ($\mu\text{g/kg}$) | As(III) ($\mu\text{g/kg}$) | As(V) ($\mu\text{g/kg}$) | iAs ($\mu\text{g/kg}$) | DMA ($\mu\text{g/kg}$) | MMA ($\mu\text{g/kg}$) |
|---------|----------------------------------|---------------------------------|-------------------------------|-----------------------------|-----------------------------|-----------------------------|
| 1568a | 290 \pm 30* | 60 \pm 12** | 39 \pm 8** | 100 \pm 20** | 171 \pm 34** | 11 \pm 2** |
| 1568 | 410 \pm 50* | 85 \pm 17** | 31 \pm 6** | 116 \pm 23** | 285 \pm 57** | 22 \pm 4** |

* Certified Value with Uncertainty expressed as a 95% Confidence Interval or 95% Confidence Interval plus an allowance for systematic error.

** Uncertainty expressed as \pm 20% of the average value from the best available data.

Standard reference materials (SRMs)

NIST (National Institute of Standards and Technology) Standard Reference Material (SRM) 1568a and 1568 rice flour were used as quality control materials for both As speciation measurements and total As determinations (Table 1). The certified values for NIST 1568, for total arsenic are $410 \pm 30 \mu\text{g}/\text{kg}$ and the certified values for As(III), DMA, and As(V) in NIST 1568 are 85 ± 17 , 285 ± 57 , and $31 \pm 6 \mu\text{g}/\text{kg}$ respectively.

Instruments

Speciation analysis was done using an Agilent 1200 Infinity LC system consisting of a vacuum degassing unit, binary pump, autosampler, and a 10-port valve (two ports were plugged). An external pump was used (LC-20AD, Shimadzu) for the internal standard spike introduction. The HPLC system was connected to the Agilent 7700x ICP-MS using $1/16 \times .0025 \times 12$ inch Peek tubing (Upchurch Scientific) and equipped with a MicroMist nebulizer and Scott double pass spray chamber. The instrument run conditions are listed in Table 2.

Standard solutions

All standard solutions were prepared daily through serial dilution of stock species standards according to the FDA method (EAM 4.11). All arsenic stock standard solutions (As(III), DMA, MMA, As(V)) were prepared in water at $1,000 \mu\text{g}/\text{mL}$. The arsenic concentrations of the DMA and MMA standards were verified using ICP-MS analysis. The total arsenic concentrations of $1 \mu\text{g}/\text{g}$ MMA and DMA standards were determined using a calibration curve of ICP-MS Arsenic Standard. These concentrations were used to recalculate the stock standard concentrations and the new values were then applied to the future calculations.

Working standards were then diluted with mobile phase to desired concentrations using the $200 \mu\text{g}/\text{mL}$ stock solution. All standard solutions were stored in the dark at 4°C . On the morning of the experiment, a fresh mixed standard solution was made. The internal standard for the speciation method was $2 \text{ ng}/\text{g}$ of As(V) in mobile phase, and it was injected post-column to monitor and to correct for signal drift. The internal standard for the totals method was rhodium/germanium with 18% isopropanol (ICP-MS Internal Standard Mix $100 \text{ mg}/\text{L}$, p/n 5188-6525).

Table 2. HPLC and ICP-MS Run Conditions

| HPLC conditions | |
|---------------------------|---|
| Guard column | Hamilton (p/n 79446) |
| Column | Hamilton PRP-X 100 anion exchange, 4.1×250 , stainless steel, $10 \mu\text{m}$ |
| Mobile phase | 10 mM ammonium phosphate pH 8.25, isocratic run |
| Flow rate | 1 mL/min |
| Injection volume | 100 μL |
| ICP-MS conditions | |
| RF power | 1,500 W |
| Plasma gas flow | 15 L/min |
| Dilution gas flow | 0.17 L/min |
| Nebulizer gas flow | 1.1 L/min |
| Nebulizer type | MicroMist |
| Sampling depth | 8.5 mm |
| Spray chamber temperature | 2°C |
| Collision cell gas | He, 4.3 mL/min |
| Data acquisition mode | Time-resolved analysis: m/z 75 for $^{75}\text{As}^+$ and m/z 77 for $^{40}\text{Ar}^{35}\text{Cl}^+$ |
| Dwell times | 0.8 s for m/z 75; 0.2 s for m/z 77 |

The standard for LOD and LOQ determination was prepared at 0.2 ng/g each of As(III), DMA, MMA, and As(V) in mobile phase. A resolution check standard was prepared using 5 ng/g As(III) and Arsenobetaine (AsB) in mobile phase, and it was used to determine the separation of AsB and As(III). A new resolution check standard was prepared when approximately 50% of As(III) had been converted to As(V).

Samples

Baby rice cereals were purchased at different locations and various grocery stores in the following US states; Illinois, California, Texas, and North Dakota. The 31 infant rice cereals were produced by seven different manufacturers. Table 3 lists all the information that could be obtained from the baby rice cereals containers.

Table 3. Summary of the Infant Rice Cereals Tested

| Brand | Sample ID | Product description | Ingredients | Location of purchase |
|-------|-----------|---|--|----------------------|
| A | Baby_1 | Mixed grain cereal | Whole wheat flour, rice whole grain oat flour, soybean oil | IL |
| A | Baby_2 | Organic brown rice | Organic whole grain brown rice flour, soy lecithin | IL |
| B | Baby_3 | Organic whole grain rice cereal | Organic whole grain brown rice, tocopherols (vitamin E), electrolytic iron | IL |
| C | Baby_4 | Organic rice cereal | Organic whole grain brown rice flour, sunflower lecithin, electrolytic iron | IL |
| A | Baby_5 | Rice single grain cereal | Organic whole grain brown rice flour, sunflower lecithin, electrolytic iron | IL |
| C | Baby_6 | Organic rice cereal | Organic whole grain brown rice flour, sunflower lecithin, electrolytic iron | TX |
| D | Baby_7 | Rice cereal, single grain | Rice flour, soybean oil, soy lecithin | CA |
| E | Baby_8 | Rice cereal | Rice flour, soybean oil, soy lecithin | ND |
| A | Baby_9 | Organic brown rice cereal | Organic whole grain brown rice, tocopherols(vit E), electrolytic iron | ND |
| B | Baby_10 | Organic brown rice cereal | Organic whole grain brown rice, tocopherols (vit E), electrolytic iron | ND |
| B | Baby_11 | Organic brown rice cereal | Organic whole grain brown rice, tocopherols (vit E), electrolytic iron | CA |
| A | Baby_12 | Organic brown rice cereal | Organic brown rice cereal | TX |
| F | Baby_13 | Rice cereal | Rice flour, sunflower lecithin | ND |
| F | Baby_14 | Rice cereal | Rice flour, sunflower lecithin | CA |
| F | Baby_15 | Rice cereal | Rice flour, sunflower lecithin | CA |
| A | Baby_16 | Rice cereal | Rice flour | TX |
| F | Baby_17 | Rice cereal | Rice flour, sunflower lecithin | TX |
| A | Baby_18 | Mixed grain cereal | Whole wheat flour, rice whole grain, oat flour, soybean oil | CA |
| A | Baby_19 | Mixed grain cereal | Whole wheat flour, rice whole grain, oat flour, soybean oil | TX |
| A | Baby_20 | Mixed grain cereal | Whole wheat flour, rice whole grain, oat flour, soybean oil | ND |
| A | Baby_21 | Mixed grain cereal | Whole wheat flour, rice whole grain, oat flour, soybean oil | CA |
| A | Baby_22 | Mixed grain cereal | Whole wheat flour, rice whole grain, oat flour, soybean oil | ND |
| A | Baby_23 | Mixed grain cereal | Whole wheat flour, rice whole grain, oat flour, soybean oil | CA |
| G | Baby_24 | Organic brown rice cereal | Organic whole grain brown rice flour, sunflower lecithin | CA |
| G | Baby_25 | Organic brown rice cereal | Organic whole grain brown rice flour, sunflower lecithin | CA |
| A | Baby_26 | Rice single grain cereal | Rice flour, tri- and dicalcium phosphate, soybean oil, soy lecitin, tocopherols, vitamins A and B, B group | CA |
| A | Baby_27 | Mixed grain cereal | Whole wheat flour, whole grain oat flour, tri- and dicalcium phosphate, soybean oil, soy lecitin, tocopherols, vitamins A and B, B group | CA |
| B | Baby_28 | Organic whole grain rice cereal with apples | Organic whole grain brown rice flour, organic apple puree, vitamin E, iron, and others | CA |
| B | Baby_29 | Organic whole grain rice cereal | Organic whole grain brown rice, vitamin E, iron | CA |
| B | Baby_30 | Organic whole grain rice cereal | Organic whole grain brown rice, vitamin E, iron | CA |
| A | Baby_31 | DHA and probiotic single grain rice cereal | Rice flour, vitamins, minerals, soybean oil, soy lecitin, tocopherols, tuna oil | CA |

Sample preparation

All baby rice cereals were analyzed for total arsenic and arsenic speciation. Infant cereals were analyzed directly from the container and were not further dried or homogenized. For the total arsenic analysis, the same hot block digestion method was used, but the samples were then injected directly into the ICP-MS without HPLC separation.

The hot block digestion method [3] was used for both the total and speciation analysis. A gravimetric method was applied to all preparations in this study. One gram of infant rice cereal sample was weighed precisely into a pre-weighed 50 mL polypropylene centrifuge tube (with lid). To this tube was added 10 g of 0.28 M HNO_3 , and the sample was vortexed for 10–30 seconds. The tightly capped tube was then placed in a preheated block digestion system at 95 °C for 90 minutes.

After the sample was cooled, 6.7 g of deionized water was added, and the rice cereal suspension was centrifuged at 3500 rpm for 10 minutes. The supernatant was passed

through a 0.45 μm nylon syringe filter attached to a 3 mL disposable syringe. Approximately the first 1 mL was discarded through the filter to waste, and 1 g of filtrate was then transferred to a 15 mL tared centrifuge tube. To this was added 2 g of pH Adjustment Solution which was prepared by adding ammonium hydroxide to the mobile phase to a pH of 9.95 ± 0.05 .

The pH of the final preparation of the method blank and samples was adjusted to between 6 and 8.5 with the pH Adjustment Solution. A portion of the resulting sample solution was transferred to a polypropylene autosampler vial for analysis by HPLC-ICP-MS

Results and Discussion

Chromatographic Separation

Figure 1 shows a typical separation. Blank tests were performed to investigate possible arsenic contamination; none was detected.




Figure 1. Chromatographic separation of arsenic standards at 5 ng/g.

Method precision

To validate the method for arsenic speciation and total methods, the NIST SRM 1568 rice flour was analyzed. Results for speciated and total As in the certified reference material were in good agreement with certified values. The repeatability of the EAM 4.11 method for analysis of baby rice cereals was determined using a set of SRM replicates (Table 4).

Repeatability (%RSD) was 1.74% for total arsenic by the As speciation method, 5.3% for the total arsenic by As total method, and 2.8% for iAs.

Calibration curves

The calibration curves for arsenic species were determined using As(III), DMA, MMA, and As(V) standards with concentrations of 0.25, 0.4, 1.0, 5.0 and 20 ng/g (Figure 2). The standard curves were linear for this range with a coefficient of variation (R^2) of 1.00. Recovery for the calibration check standard, which was composed of 10 ng/g each of As(III), DMA, MMA, and As(V), was 99–113% with % RSD < 4%.

Table 4. Method Precision

| Sample name (separate preparations) | As(III) ($\mu\text{g}/\text{kg}$) | DMA ($\mu\text{g}/\text{kg}$) | MMA ($\mu\text{g}/\text{kg}$) | As(V) ($\mu\text{g}/\text{kg}$) | iAs ^a ($\mu\text{g}/\text{kg}$) | Total As (by speciation method) | Total As (by totals method) | Mass balance ^b |
|--|--|------------------------------------|------------------------------------|--------------------------------------|---|------------------------------------|--------------------------------|------------------------------|
| Certified value | 85 \pm 17 | 285 \pm 57 | 22 \pm 4 | 31 \pm 6 | 116 \pm 23 | 410 \pm 50 | 410 \pm 50 | |
| 1568_rep1 | 65 | 289.9 | 24.8 | 40.6 | 105.6 | 420.3 | 398.4 | 105.5 |
| 1568_rep2 | 68.5 | 289 | 24.6 | 42.5 | 110.9 | 424.5 | 445.1 | 95.4 |
| 1568_rep3 | 66.9 | 300.6 | 25 | 44.7 | 111.6 | 437.2 | 405.2 | 107.9 |
| 1568_rep4 | 70.3 | 288.7 | 21.6 | 37.5 | 107.8 | 418.1 | 442 | 94.6 |
| 1568_rep5 | 69.8 | 289.4 | 22.8 | 43.3 | 113.1 | 425.4 | 442.3 | 96.2 |
| Average | 68.1 | 291.5 | 23.8 | 41.7 | 109.8 | 425.1 | 426.6 | 99.9 |
| Standard deviation (SD) | 2.2 | 5.1 | 1.5 | 2.8 | 3.1 | 7.4 | 22.8 | 6.3 |
| %RSD | 3.2 | 1.7 | 6.3 | 6.7 | 2.8 | 1.7 | 5.3 | 6.3 |

^a iAs equals the sum of As(III) and As(V).

^b Mass balance = Total As from speciation (As(III) + DMA + MMA + As(V)) / (total As from total As method)




Figure 2. Calibration curves for the EAM 4:11 method.

Limits of detection and quantitation

The Analytical Solution Detection Limit (ASDL), Analytical Solution Quantitation Limit (ASQL), and method Limit of Detection (LOD) and Limit of Quantitation (LOQ) were determined for the arsenic speciation analysis in rice using the HPLC-ICP-MS method (Table 5). The LOD values for As species ranged from approximately 0.9 ng/g to 1.8 ng/g, and LOQ from 7 µg/kg to 14 µg/kg.

Recoveries

Generally, spiked sample recoveries for iAs, DMA and MMA were excellent, considering arsenic conversion from As(III) to As (V). The recoveries for iAs, DMA, and MMA ranged from a low of 97.3% for iAs to a high of 115.0% for DMA (Table 6).

Table 6. % Recovery of Spiked Arsenic (n = 2)

| Spike concentration of each species | As(III) | DMA | MMA | As(V) | iAs |
|-------------------------------------|--------------|-------------|-------------|--------------|-------------|
| 50 ng/g | 104.1 ± 38.2 | 115.0 ± 1.2 | 104.8 ± 1.1 | 109.1 ± 38.3 | 106.5 ± 0.6 |
| 75 ng/g | 85.1 ± 1.3 | 108.6 ± 2.9 | 106.8 ± 0.8 | 116.4 ± 0.9 | 100.5 ± 1.1 |
| 100 ng/g | 67.8 ± 32.8 | 112.4 ± 2.2 | 105.0 ± 0.8 | 143.7 ± 28.5 | 105.2 ± 2.6 |
| 150 ng/g | 74.7 ± 10.1 | 104.0 ± 9.4 | 106.0 ± 6.1 | 120.7 ± 7.7 | 97.3 ± 8.9 |
| 150 ng/g | 67.9 ± 15.6 | 112.9 ± 2.8 | 106.6 ± 0.5 | 142.8 ± 12.4 | 104.8 ± 1.8 |
| 225 ng/g | 88.1 ± 1.5 | 109.3 ± 1.2 | 105.2 ± 1.0 | 118.2 ± 3.0 | 102.9 ± 2.2 |

Total and speciation arsenic concentrations in baby rice cereals

Table 7 shows the total and speciation results for arsenic in baby rice cereals. There was no significant difference between the two duplicates run on each cereal sample. All 31 of the cereals tested contained arsenic, as well as inorganic arsenic (iAs), which is defined as the sum of the As(II) and As(V) concentrations. The average total arsenic and iAs concentrations were 174.4 and 101.4 µg/kg, respectively in rice cereals. The average percentage of iAs as a component of the average total arsenic in the baby rice cereals tested was 58.2%.

The mass balance between the As speciation and totals methods (sum of all As species by the speciation method, divided by the As concentration from the totals method) ranged from 85.8% to 106.0%, with an average mass balance of 96.7%. All mass balances are shown in Table 7.

Table 5. ASDL, ASQL, LOD and LOQ of Arsenic Speciation Analysis by HPLC-ICP-MS

| | As(III) | DMA | MMA | As (V) |
|-------------|-------------|-------------|-------------|-------------|
| ASDL (ng/g) | 0.034 | 0.017 | 0.024 | 0.035 |
| LOD (µ/kg) | 1.68 | 0.87 | 1.18 | 1.74 |
| ASQL (ng/g) | 0.262 | 0.135 | 0.184 | 0.273 |
| LOQ (µ/kg) | 13.1 | 6.8 | 9.2 | 13.6 |

The average total arsenic and iAs concentrations were lower in mixed grain cereals (105.1 and 62.8 µg/kg, respectively). The average percentage of iAs as a component of the average total arsenic of the average total arsenic in the baby mixed grain cereals tested was thus 59.8%.

As listed in Table 3, all 31 baby rice cereals were made from the flour of either organic rice, organic rice whole grain, mixed grain (whole wheat, rice whole grain, and oat), or rice. Results show that levels of inorganic arsenic (iAs) varied greatly across all rice cereals (Table 7). There was no substantial difference in iAs levels between organic and conventional rice cereals (an average of 98 ng/g in organic white rice cereal versus 101 ng/g in conventional white rice cereal).

In the baby rice cereals tested, the total arsenic concentration ranged from 79.9 to 277 ng/g. The iAs, as a percent of total arsenic, ranged from 23 to 81%. The major organic arsenic species detected was DMA (Figure 3). MMA was either not detected or only trace levels were found (Table 7).

Table 7. Results from the 23 Infant Rice Cereals Analyzed in this Application Note

| Sample ID | As(III) (µg/kg) | DMA (µg/kg) | MMA (µg/kg) | As(V) (µg/kg) | iAs ^a (µg/kg) | Total As (by speciation) (µg/kg) | Total As (by totals) (µg/kg) | % Mass balance ^b | Inorganic As per serving (µg) |
|--|--------------------|----------------|----------------|------------------|-----------------------------|--|------------------------------------|--------------------------------|-------------------------------------|
| Mixed grain cereal_1 | 20.6 ± 7.7 | 27.5 ± 0.2 | TR (1.4) | 36.6 ± 8.4 | 57.2 ± 0.7 | 86.2 ± 0.5 | 88.3 ± 1.3 | 97.6 | 0.9 |
| Mixed grain cereal_18 | 18.1 ± 3.3 | 35.2 ± 1.4 | TR (1.9) | 44.9 ± 4.1 | 63.0 ± 0.7 | 100. ± 2.0 | 112. ± 5.7 | 89.1 | 0.9 |
| Mixed grain cereal_19 | 12.1 ± 9.2 | 23.0 ± 0.7 | TR (1.5) | 43.3 ± 7.9 | 55.5 ± 1.3 | 79.9 ± 1.4 | 88.0 ± 1.4 | 90.9 | 0.8 |
| Mixed grain cereal_20 | 4.84 ± 1.3 | 50.5 ± 0.6 | TR (2.3) | 58.9 ± 1.1 | 63.7 ± 0.1 | 116. ± 0.6 | 124. ± 0.5 | 93.8 | 1.0 |
| Mixed grain cereal_21 | 13.0 ± 2.1 | 48.2 ± 0.4 | TR (2.4) | 58.5 ± 2.0 | 71.6 ± 0.1 | 122. ± 0.5 | 126. ± 0.4 | 96.4 | 1.1 |
| Mixed grain cereal_22 | 4.83 ± 0.5 | 50.0 ± 1.7 | TR (2.5) | 64.0 ± 1.1 | 68.8 ± 0.6 | 121. ± 2.7 | 131. ± 9.1 | 92.2 | 1.0 |
| Mixed grain cereal_23 | 9.89 ± 9.6 | 49.1 ± 0.6 | TR (3.0) | 55.4 ± 9.1 | 65.3 ± 0.4 | 117. ± 0.7 | 128. ± 4.6 | 91.5 | 1.0 |
| Mixed grain cereal_27 | 19.3 ± 0.2 | 38.1 ± 0.5 | TR (1.3) | 37.7 ± 5.3 | 57.1 ± 5.0 | 96.5 ± 3.9 | 109. ± 0.4 | 88.1 | 0.9 |
| Organic, rice cereal_4 | 37.8 ± 8.6 | 153. ± 0.6 | TR (8.9) | 67.3 ± 4.9 | 105. ± 3.7 | 267. ± 3.8 | 272. ± 7.9 | 98.2 | 1.6 |
| Organic, rice cereal_6 | 26.0 ± 0.4 | 51.8 ± 1.5 | TR (2.8) | 64.7 ± 1.9 | 90.7 ± 1.4 | 145. ± 3.1 | 146. ± 7.3 | 99.3 | 1.4 |
| Rice single grain cereal_5 | 65.7 ± 5.7 | 100. ± 0.1 | TR (4.9) | 72.0 ± 3.7 | 137. ± 1.9 | 243. ± 2.5 | 230. ± 10. | 106 | 2.1 |
| Rice single grain cereal_7 | 32.9 ± 0.8 | 89.2 ± 0.9 | TR (5.1) | 85.9 ± 0.3 | 118. ± 1.1 | 213. ± 0.7 | 212. ± 8.5 | 100 | 1.8 |
| Rice single grain cereal_8 | 56.6 ± 1.8 | 126. ± 1.0 | TR (6.6) | 76.1 ± 0.5 | 132. ± 1.3 | 265. ± 0.8 | 256. ± 0.4 | 104 | 2.0 |
| Rice single grain cereal_13 | 27.6 ± 0.1 | 35.3 ± 0.0 | TR (1.5) | 45.2 ± 1.7 | 72.8 ± 1.6 | 109. ± 1.2 | 108. ± 1.7 | 101 | 1.1 |
| Rice single grain cereal_14 | 24.7 ± 1.0 | 32.2 ± 0.9 | TR (3.0) | 43.5 ± 0.0 | 68.2 ± 0.9 | 103. ± 1.8 | 101. ± 1.7 | 102 | 1.0 |
| Rice single grain cereal_15 | 24.1 ± 1.6 | 30.0 ± 0.3 | TR (1.6) | 33.8 ± 0.0 | 58.0 ± 1.7 | 89.7 ± 3.1 | 92.5 ± 3.2 | 97.1 | 0.9 |
| Rice single grain cereal_16 | 59.9 ± 22. | 107. ± 1.5 | TR (3.4) | 65.6 ± 14. | 125. ± 8.1 | 236. ± 9.6 | 227. ± 1.6 | 104 | 1.9 |
| Rice single grain cereal_17 | 35.2 ± 10. | 34.3 ± 0.9 | TR (1.8) | 54.7 ± 7.0 | 90.0 ± 3.1 | 126. ± 3.1 | 127. ± 1.4 | 99.2 | 1.4 |
| Rice single grain cereal_26 | 60.1 ± 8.3 | 97.1 ± 1.8 | TR (2.7) | 47.4 ± 8.2 | 107. ± 0.1 | 207. ± 2.1 | 241. ± 2.4 | 85.8 | 1.6 |
| Rice single grain cereal_31 | 58.6 ± 3.7 | 104. ± 1.6 | TR (2.8) | 39.1 ± 2.7 | 97.8 ± 0.9 | 205. ± 3.5 | 237. ± 0.9 | 86.4 | 1.5 |
| Organic, whole grain rice cereal_3 | 16.2 ± 7.3 | 15.9 ± 0.5 | 0.0 | 51.7 ± 5.4 | 68.0 ± 1.9 | 84.4 ± 3.0 | 85.1 ± 1.7 | 99.1 | 1.0 |
| Organic brown rice cereal_2 | 21.7 ± 0.2 | 42.9 ± 0.1 | TR (4.2) | 111. ± 1.2 | 133. ± 1.0 | 180. ± 0.2 | 178. ± 0.1 | 101 | 2.0 |
| Organic brown rice cereal_9 | 21.7 ± 10. | 27.7 ± 1.0 | 0.0 | 34.3 ± 9.1 | 56.0 ± 1.7 | 84.7 ± 2.6 | 83.9 ± 2.6 | 101 | 0.8 |
| Organic brown rice cereal_10 | 24.2 ± 1.3 | 206. ± 4.2 | TR (6.5) | 39.5 ± 2.6 | 63.8 ± 4.0 | 277. ± 8.0 | 261. ± 7.0 | 106 | 1.0 |
| Organic brown rice cereal_11 | 21.3 ± 3.1 | 38.5 ± 0.3 | TR (2.8) | 107. ± 0.3 | 128. ± 2.7 | 169. ± 2.0 | 165. ± 7.8 | 103 | 1.9 |
| Organic brown rice cereal_12 | 23.0 ± 3.5 | 107. ± 1.0 | TR (3.4) | 65.7 ± 1.3 | 88.8 ± 2.2 | 200. ± 2.9 | 206. ± 13. | 97.2 | 1.3 |
| Organic brown rice cereal_24 | 18.0 ± 0.9 | 25.8 ± 0.0 | TR (1.7) | 69.4 ± 0.7 | 87.4 ± 1.6 | 114. ± 2.0 | 127. ± 1.9 | 90.2 | 1.3 |
| Organic brown rice cereal_25 | 23.1 ± 1.8 | 27.6 ± 0.0 | TR (2.0) | 73.0 ± 0.6 | 96.2 ± 2.4 | 125. ± 1.9 | 133. ± 0.0 | 94.2 | 1.4 |
| Organic, whole grain rice cereal with apples_28 | 41.8 ± 1.2 | 32.2 ± 0.0 | TR (1.5) | 63.3 ± 0.4 | 105. ± 0.7 | 139. ± 0.9 | 160. ± 1.0 | 86.6 | 1.6 |
| Organic, whole grain rice cereal_29 | 51.5 ± 2.2 | 65.0 ± 0.3 | TR (2.6) | 89.6 ± 1.2 | 141. ± 0.9 | 208. ± 0.1 | 242. ± 4.4 | 86.0 | 2.1 |
| Organic, whole grain rice cereal_30 | 81.7 ± 7.8 | 51.8 ± 0.6 | TR (1.2) | 76.7 ± 13. | 158. ± 6.0 | 211. ± 5.4 | 244. ± 0.0 | 86.4 | 2.4 |

^a iAs equals the sum of As(III) and As(V)

^b Mass balance: Total As from speciation (As(III) + DMA + MMA + As(V)) / total As from totals method; all samples were analyzed in duplicate (n = 2).

TR = trace level




Figure 3. A representative chromatogram overlay of arsenic speciation in five infant rice cereals.

The mixed grain rice cereal contained the least total (105 ng/g average value) and inorganic arsenic (63 ng/g average value), by speciation.

The results of the analysis of the baby rice cereals were compared with results obtained in the literature [1,8]. Concentrations of total arsenic and inorganic arsenic were not as high as have been reported for some of the infant and toddler formulas. This is due to the fact that none of the baby rice cereals in this study contained organic brown rice syrup (OBRS). It has been reported that toddler formulas containing OBRS as the primary ingredient can contain arsenic levels over 20 times higher than those in non-OBRS formulas [1].

Comparison to FDA findings

On September 19, 2012, the FDA released the first analytical results of nearly 200 samples of rice and rice products collected in the U.S. marketplace. Among these samples, 32 were rice cereal samples and only 3 were infant rice cereal samples (Table 8). The FDA reported iAs values were much higher than our results. Average inorganic arsenic (iAs) was 3.5 µg/serving in FDA tested rice cereals (range of 1.5–9.7 µg/serving), and 2.7 µg/serving in infant rice cereals.

The 23 infant rice cereal samples that were analyzed in this study had an average 1.5 µg of iAs per serving, and eight mixed grain cereals contained 0.9 µg/serving iAs on average. The iAs values ranged from 0.8 to 2.4 µg/serving (Table 7).

Table 8. FDA Summary Analytical Results for Infant Rice Cereals, from the "Rice/Rice Product Sampling - September 2012"

| Sample ID | Product category | Sample description | Country of origin | Inorganic As per serving (µg/serving*) |
|-----------|------------------|--|-------------------|--|
| 70145 | Rice Cereal | Organic whole grain rice cereal (infant) | NR | 3.2 |
| 720341 | Rice Cereal | Organic whole grain rice cereal (infant) | NR | 2.9 |
| 719983 | Rice Cereal | Rice single grain (infant) | NR | 2.0 |

*Serving size based on Reference Amount Customarily Consumed (RACC) per 21CFR 101.12 for infant rice cereals, which used a value of 15 g to calculate the µg of inorganic arsenic per serving.

Conclusions

This study demonstrated the successful application of the HPLC-ICP-MS method specified in the FDA Elemental Analysis Manual (EAM): Section 4.11. for the determination of arsenic species in infant rice cereals. The sensitivity (LODs and LOQs) and precision (%RSD) were sufficient to detect low levels of all arsenic species, including iAs levels that were below those previously found in baby rice cereals in an FDA study. This study provides new and much needed information for the assessment of arsenic levels in rice and rice-based infant cereals.

Acknowledgments

The authors wish to acknowledge Amir Liba and Mike Woodman of Agilent Technologies for technical support and Agilent Technologies for provision of the HPLC-ICP-MS equipment.

For More Information

These data represent typical results. For more information on our products and services, visit our Web site at www.agilent.com/chem.

References

1. B. Jackson, *et al.*, "Arsenic, Organic Foods, and Brown Rice Syrup", *Environ. Health Perspect.*, 2012, 120, 623-626.
2. B. Jackson, *et al.*, "Arsenic concentration and speciation in infant formulas and first foods", *Pure Appl. Chem.*, 2012, 84(2), 215-225.
3. A. Carey, *et al.*, "Phloem transport of arsenic species from flag leaf to grain during grain filling", *New Phytol.* 2011, 192(1), 87-98.
4. A. Meharg, *et al.*, "Arsenic & Rice", 2012, ISBN-13, 978-9400729469
5. A. Meharg, *et al.*, "Geographical Variation in Total and Inorganic Arsenic Content of Polished (White) Rice", *Environ. Sci. Technol.*, 2009, 43(5), 1612-1617.
6. G. Sun, *et al.*, "Inorganic arsenic in rice bran and its products are an order of magnitude higher than in bulk grain", *Environ. Sci. Technol.*, 2008, 42(19), 7542-7546.
7. R. Juskelis, W. Li, J. Nelson, J.C. Cappozzo, "Arsenic speciation in rice cereals for infants", *J. Agric. Food Chem.*, 2013, DOI: 10.1021/jf401873z
8. N. Vela, D. Heitkemper, "Total arsenic determination and speciation in infant food products by ion chromatography- inductively coupled plasma- mass spectrometry", *J. AOAC Int.*, 2004, 87(1), 244-252.

www.agilent.com/chem

Agilent shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

Information, descriptions, and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc., 2013
Printed in the USA
November 8, 2013
5991-2568EN



Agilent Technologies