Azura

▶ **Autosampler 3950
User Manual**

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V6818

HPLC

Table of Contents

Note For your own safety, **read** the manual and **always** observe the warnings and safety information on the device and in the manual!

Intended Use	7
Device Overview	7
Features	7
Standard	7
Device Variants	8
Solvents	9
Flushing Solution	10
Safety for Users	11
Definition of Personal and Material Damage	13
Decontamination	13
Symbols and Signs	14
Installation	15
Scope of Delivery	15
Checking the Scope of Delivery	15
Unpacking and Setup	16
Packaging and Transport	16
Contacting the Technical Support	16
Location Requirements	16
Space requirements	17
Unpacking	17
System Adapter	18
Start-Up	19
Front View of the Device	19
Front View with Sample Compartment	21
Rear View of the Device	22
Local Area Network and Automatic Configuration	23
Connecting the Device in a Local Area Network (LAN) to a Computer	23
Injection Principles	25
ILD™ of Spark Holland	25
PASA™ Loop Injection Principle	25
Injection Methods	25
Full Loop Injection	27
Air Segment for Full Loop Injection	29
Partial Loop Fill	30
Microliter Pick-Up Injection	33
Air Segment with µl-Pick-Up Injection	35
Microliter Pick-Up Injection with 84+3	36
Air Segment with µl-Pick-Up Injection 84+3	37
84+3 Vial Plate	39

<ul-pick-up .="" .<="" 84+3="" for="" parameters="" plate="" td="" the="" vial=""> <td>40</td> </ul-pick-up>	40
Details on Programming the 84+3 Injection Method	40
Air Needles	42
Standard Air Needle	42
Choosing the Correct Air Needle	43
Calculation Example for Air Needle	45
Handling the Sample Vials	45
Mixing and Diluting	46
An example: Add	46
An example: Mix	47
Sample Positions in Mix Method	47
Processing in Columns	47
Processing in Rows	48
Parameters for Mixing Method with 84+3 Vial Plate	48
Details on Programming the 84+3 Mixing Method	49
Capillary and Hose Connections	50
Connecting the Injection Valve	51
Connecting the Syringe	51
Tube Guide for Flushing Solution	51
Connecting the Drainage Tubing	53
PEEK Connections	53
Connecting the Autosampler to Other Devices	54
Controlling the Autosampler with Chromatography Software	54
Checking and Configuring the Parameters of the Autosampler	54
Configuration Window of ClarityChrom®	54
Autosampler Device Software	55
Flushing the System	55
System Flushing with Autosampler 3950 Service Manager	55
I/O Connection	56
Defining the TTL Inputs	56
Defining the Closed-Contact Output	57
Configuration of I/O Connection (9 pins)	57

Device Test	58
Test Intervals	58
Devices and Components for the Test	58
1. Reproducibility of Sample Volume	59
Standard Setting of Autosampler	59
Method Parameters of Pump	59
Method Parameters of Autosampler	59
Method Parameters of UV Detector	59
Configuring Repeat Runs of Autosampler	59
Starting Repeat Runs	59
Analyzing the Individual Chromatograms	60
Formula for Determining the Arithmetic Mean	60
2. Sample Carryover	60
Creating a Sequence with 6 Lines	60
Analyzing the Individual Chromatograms	61
Formula for Calculating Sample Carryover	61
3. Linearity	61
Analyzing the Individual Chromatograms	61
Formula for Determining the Correlation	62
4. Mixture Test	62
Creating a Sequence with 2 Lines	62
Positioning Vials for Dilution	63
Analyzing the Individual Chromatograms	63
Archiving	63
Test Report	64
Operation Qualification Report	65
Maintenance and Care	66
Contacting the Technical Support	66
Maintenance Contract	66
Which Type of Maintenance Tasks May Users Perform on the Device?	66
Leaks in the Capillary Screw Fittings	66
Exchanging the Fuses	66
Changing the Valve and Rotor Seal	67
Removing the Injection Valve and Rotor Seal	68
Installing the Injection Valve	69
System Flushing	69
Exchanging the Sample Loop	70
Exchanging the Sample Needle	71
Exchanging the Air Needle	72
Exchanging the Syringe	73
Exchanging the Syringe Plunger or Plunger Tip	74
Exchanging the Syringe Valve	75
Cleaning and Caring for the Device	76
Disposal	77

Troubleshooting	78
Device Errors	78
Checking the Valve	78
LAN	78
Analytical Errors	79
System Messages from OpenLAB®	82
Technical Data	94
Standard Version	94
Product Range	95
Device and Accessories	95
Repeat Orders	95
Accessories 84 +3	96
Legal Information	97
Warranty Conditions	97
Transport Damage	97
HPLC Glossary	98
List of Figures	99
Index	101

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- support@knauer.net
- +49 30 8015010

Vielen Dank.

To whom it may concern In case you prefer a French language user manual for this product, submit your request including the corresponding serial number via email or fax to KNAUER:

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- support@knauer.net
- +49 30 8015010

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Intended Use

Note Only use the device for applications that fall within the range of the intended use. Otherwise, the protective and safety equipment of the device could fail.

Device Overview

Autosampler 3950 The Autosampler 3950 was developed for chromatographic analyses in the high pressure (HPLC) range. An injection valve guarantees precise injection volumes with a maximum of 1000 bar. Additional variants offer temperature control. The device features fast sample injection, rapid flushing cycles and high sample throughput. The Autosampler is generally installed as a core element in an HPLC system.

Legend

- ① Opening for capillary feed
- ② Removable cover
- ③ Inward sliding glass door
- ④ Temperature control with coolunit inside the Autosampler
- ⑤ Hose connectors


Fig. 1 Autosampler 3950, with optional temperature control and coolunit

Scope of application The device can be used in the following areas:

- Biochemistry analysis
- Chemical analysis
- Food analysis
- Pharmaceutical analysis
- Environmental analysis

Features

Standard

- Injection valve with injection volume maximum of 1000 bar
- Injection volume for HPLC between 0.1 µl of maximum 5000 µl

- Intermediate Loop Decompression (ILD™) technology¹
- Temperature control between 4 – 40 °C²
- Flexibly equipped with microtiter plate or standard sample plate
- Maximum sample capacity of 768 Wells or 96 standard Autosampler vials
- Injection cycle < 60 s, 60 s incl. cleaning
- Full loop or partial loop injection or microliter pick-up injection
- Quick-exchange injection valve
- Sample priority function
- High-resolution syringe controller

Device Variants

Two variants of the device are available:

- Autosampler without temperature control
- Autosampler with temperature control

1. ILD™ is a trademark of Spark Holland
2. Relevant to Autosampler 3950 with temperature control

Solvents

Even small quantities of other substances, such as additives, modifiers, or salts can influence the durability of the materials.

If there is any doubt, contact the Technical Support of the manufacturer.

Suitable solvents

- Acetate buffer solutions
- Acetone at 4 °C–25 °C (39.2 °F–77.0 °F)¹
- Acetonitrile²
- Benzene
- Carbon dioxide (liquid 99.999 % CO₂)
- Chloroform
- Dilute acetic acid (e.g. 0.1–1 %) at 25 °C/77.0 °F
- Dilute ammonia solution
- Dilute sodium hydroxide (1 M)
- Ethyl acetate
- Ethanol
- Formiate buffer solution
- Isopropanol
- Methanol
- Phosphate buffer solutions (0.5 M)
- Toluol
- Water

1. valid for the specified temperature range

2. not recommended in combination with PEEK small parts and PEEK capillary

Less suitable solvents

- Diethylamine (0.1 %) (DEA)
- Dilute phosphoric acid
- Dimethyl sulfoxide (DMSO)
- Methylene chloride¹
- Slightly volatile solvents
- Tetrahydrofuran (THF)¹
- Triethylamine (0.1 %) (TEA)
- Trifluoroacetic acid (0.1 %) (TFA)

1. not recommended in combination with PEEK small parts and PEEK capillary

Not suitable solvents

- Concentrated mineral and organic acids
- Concentrated bases
- Halogenated hydrocarbons, e.g. Freon[®]
- Perfluorinated solvents, e.g. Fluorinert[®] FC-75, FC-40
- Perfluorinated polyether, e.g. Fomblin[®]
- Solvents containing particles

Solvents not suitable for degassers

- Azides
- Benzene
- Carbon dioxide (liquid 99.999 % CO₂)
- Concentrated mineral and organic acids
- Concentrated bases
- Dilute sodium hydroxide (1 M)
- Halogenated hydrocarbons, e.g. Freon[®]
- Hexafluoroisopropanol
- Hexanes (60 % *n*-Hexane)
- Hydro fluoro solvents
- Perfluorinated solvents, e.g. Fluorinert[®] FC-75, FC-40
- Perfluorinated polyether, e.g. Fomblin[®]
- Solvents containing particles

Flushing Solution

Do not use a salt or buffer solution for flushing!

Safety for Users

- Professional group** The user manual addresses persons who are qualified as chemical laboratory technicians or have completed comparable vocational training.
- The following knowledge is required:
- Fundamental knowledge of liquid chromatography
 - Knowledge regarding substances that are suitable only to a limited extent for use in liquid chromatography
 - Knowledge regarding the health risks of chemicals
 - Participation during an installation of a device or a training by the company KNAUER or an authorized company.
- If you do not belong to this or a comparable professional group, you may not perform the work described in this user manual under any circumstances. In this case, please contact your superior.
- Safety equipment** When working with the device, take measures according to lab regulations and wear protective clothing:
- Safety glasses with side protection
 - Protective gloves
 - Lab coat
- What must be taken into account?**
- All safety instructions in the user manual
 - The environmental, installation, and connection specifications in the user manual
 - National and international regulations pertaining to laboratory work
 - Original spare parts, tools, and solvents made or recommended by KNAUER
 - Good Laboratory Practice (GLP)
 - Accident prevention regulations published by the accident insurance companies for laboratory work
 - Filtration of substances under analysis
 - Use of inline filters
 - Once they have been used, never re-use capillaries in other areas of the HPLC system.
 - Only use a given PEEK fitting for one specific port and never re-use it for other ports. Always install new PEEK fittings on each separate port.
 - Follow KNAUER or manufacturer's instructions on caring for the columns

More safety-relevant information is listed in alphabetical order in the following table:

Topic	Explanations
Flammability	Organic solvents are highly flammable. Since capillaries can detach from their screw fittings and allow solvent to escape, it is prohibited to have any open flames near the analytical system.
Solvent tray	Risk of electrical shock or short circuit if liquids get into the device's interior. For this reason, place all bottles in a solvent tray.
Solvent lines	Install capillaries and tubing in such a way that liquids cannot get into the interior in case of a leak.
Leaks	Regularly check if any system components are leaking.
Power cable	Defective power cables are not to be used to connect the device and the power supply system.
Self-ignition point	Only use eluents that have a self-ignition point higher than 150 °C under normal ambient conditions.
Power strip	If several devices are connected to one power strip, always consider the maximum power consumption of each device.
Power supply	Only connect devices to voltage sources, whose voltage equals the device's voltage.
Toxicity	Organic eluents are toxic above a certain concentration. Ensure that work areas are always well-ventilated! Wear protective gloves and safety glasses when working on the device!

Where is use of the device prohibited?

Never use the system in potentially explosive atmospheres without appropriate protective equipment. For further information, contact the Technical Support of KNAUER.

Taking the device out of operation

At any time, take the device completely out of operation by either switching off the power switch or by pulling the power plug.

Opening the device

The device may be opened by the KNAUER Technical Support or any company authorized by KNAUER only.

Definition of Personal and Material Damage

Possible dangers related to the device are divided into personal and material damage in this user manual.

Category	Explanations
DANGER!	Lethal or very serious injuries can occur.
WARNING!	Serious or moderate injuries can occur.
CAUTION!	Moderate injuries can occur. Device defects can occur.

Decontamination

Contamination of devices with toxic, infectious or radioactive substances poses a hazard for all persons during operation, repair, sale, and disposal of a device.

DANGER!

Danger if getting in contact with toxic, infectious or radio-active substances.

Before disposing off or sending away contaminated devices, commission an expert with the decontamination.

All contaminated devices must be properly decontaminated by a specialist company or the operating company before they can be recommissioned, repaired, sold, or disposed of.

All materials or fluids used for decontamination must be collected separately and disposed of properly.

Symbols and Signs

The following symbols and signs can be found on the device, in the chromatography software or in the user manual:

	Symbol	Meaning
Warning signs		Electric shock hazard
		Danger caused by potentially toxic substances
		Heavy weight hazard, back injuries can occur.
		Sharp edges hazard, stab wounds can occur.
		General warning sign, moderate injuries can occur and also damages to system, device, or components.
		Electrostatic discharge hazard, damages to system, device, or components can occur.
Mandatory signs		Wear protective gloves to avoid skin irritations.
CE mark		A device or system marked with CE fulfills the product specific requirements of European directives. This is confirmed in a Declaration of Conformity.
		Testing seals in Canada and the USA at nationally recognized testing centers (NRTL). The certified device or system has successfully passed the quality and security tests.

Installation

Scope of Delivery

Note Only use original parts and accessories made by KNAUER or a company authorized by KNAUER.

Delivery	Autosampler and 250 µl syringe	<input type="checkbox"/>
	User manual German/English	<input type="checkbox"/>
	<ul style="list-style-type: none"> ▪ Cable: <ul style="list-style-type: none"> ▪ Power cable ▪ Network cable ▪ I/O interface cable, 9 pins 	<input type="checkbox"/>
	CD <i>Autosampler 3950 Service Manager</i>	<input type="checkbox"/>
Accessories kit	Sample loop, 10 µl and 100 µl	<input type="checkbox"/>
	PTFE tube, 300 cm in length	<input type="checkbox"/>
	Silicon drainage tube, 200 cm	<input type="checkbox"/>
	2 x vial plate for 48 vials, 1.5 ml	<input type="checkbox"/>
	5 x microtiter plate, 96 Wells	<input type="checkbox"/>
	2 x 2.5 A fuses	<input type="checkbox"/>
	1 x wash bottle	<input type="checkbox"/>
	Optional: Temperature control	<input type="checkbox"/>

Checking the Scope of Delivery

1. Check whether the device and accessories are complete.
2. If a part is missing, inform the technical support of KNAUER.

Unpacking and Setup

Packaging and Transport

At the factory, the device is carefully packed for safe transport.

Note Check for damage caused during transportation. In case you notice any damage, contact the technical support and the forwarder company.

Contacting the Technical Support

You have various options to contact the Technical Support:

Phone +49 30 809727-111

Fax +49 30 8015010

E-mail support@knauer.net

You can make your requests in English and German.

Location Requirements

Requirements The location for the device must meet the following requirements:

- protect from heavy ventilation
- Weight of the Autosampler 21 kg (with temperature control)
- Dimensions (width × height × depth) 300 × 377 × 577 mm
- Line voltage 95–240 V DC
- Humidity 20–80 % RH
- Temperature 10–40 °C
50–104 °F

CAUTION!

Defect of the device due to overheating possible! Protect the device against exposure to direct sunlight.

Make sure the room is well-ventilated.

Keep clear at least 15 cm at rear and 5 – 10 cm at each side for air circulation.

Note The autosampler is designed exclusively for use in closed areas with an ambient temperature of 10–40 °C.

Space requirements

Side clearance to other devices:

- If there is a device on one side, min. clearance of 5 cm.
- If there are devices on both sides, min. clearance of 10 cm.

Unpacking

WARNING! Back injuries while lifting or carrying the device caused by its heavy weight. To avoid injuries, you should ask a second person for help.

Store all packing material. Retain included packing list carefully for repeat orders.

Tool	Utility knife
Procedure	<ol style="list-style-type: none"> 1. Setup the delivery in such a way that the label is in the correct position. Using the utility knife, cut the adhesive tape. Open the packaging. 2. Remove the foam insert. Take out the accessories kit and the manual. 3. Open the accessories kit and take out all accessories. Check the scope of delivery. In case any parts are missing, contact the Technical Support. 4. Clasp the device at its side panels and lift it out of the packaging. 5. Remove the foam inserts from the device. 6. Check for damage caused during transportation. In case you notice any damage, contact the Technical Support. 7. Set-up the device in its location. 8. Remove the adhesive tape from the door of the autosampler that was used as transportation protection.

System Adapter

The autosampler can be integrated into an HPCL system with the help of a system adapter. The system adapter is a metal sheet which can be mounted onto the top cover of the autosampler. Afterwards, additional AZURA devices can be set up on the autosampler.


Fig. 2 System adapter for autosampler

Prerequisite The screws at the autosampler have been loosened and removed.

Tool Allen wrench, 2.5 mm

CAUTION! Too much weight can cause damage to the device! Do not load the system adapter with more than 60 kg!

Process and figure

1. Remove the screws from the autosampler left and right of the front ②.
2. Place the system adapter onto the autosampler. The bore holes point to the front of the autosampler ①.
3. Attach the system adapter by tightening the screws ②.


Fig. 3 System adapter installed

Next steps Place the devices on the system adapter.

Start-Up

Note Before initial startup, wait approximately one hour until the temperature of the device has adapted to the ambient temperature.


Auto-Injection System The speed of the auto-injection system has been increased to fulfill the requirements of ultra-high performance liquid chromatography. However, be aware that the high speed of the auto-injection system can cause stab injuries when handled inappropriately!

Front View of the Device

The door of the Autosampler can be pushed horizontally into the interior of the device. The front paneling can be completely removed. In order to replace the vial plate, the cooler cover has to be removed.

CAUTION! Personal injuries can occur when the front panel is open or has been removed. Always operate the autosampler with its front panel closed.

Opening the door Position your hands on both sides of the door ① and open. As shown in the figure, push the door inwards.


Fig. 5 Push the door into the interior

Removing the front panel Press the buttons on the side ① at the same time and remove the front paneling.


Fig. 6 Removing the front panel

Removing the coolunit cover Pull the coolunit cover in the direction indicated by the arrow.



Fig. 7 Removing the coolunit cover

Front View with Sample Compartment

Legend

- ① Valve
- ② Syringe
- ③ Needle guide
- ④ Collecting container
- ⑤ Flushing solution bottle
- ⑥ Tube connector for waste liquid
- ⑦ Tube connector for condensed water

Rear View of the Device

The following components, connections and warnings can be found at the rear panel of the device:

Legend

- ① LAN port
- ② I/O connection (9-pin)
- ③ Power switch
- ④ CE mark
- ⑤ Fuse box
- ⑥ Mains power connection
- ⑦ Cooling fan (optional)
- ⑧ Serial number and year of manufacture of device


Fig. 9 Rear panel of the device

Local Area Network and Automatic Configuration

The autosampler is operated exclusively using the chromatography software.

- Remote control** Normally, the autosampler is controlled by the chromatography software through a local network (LAN).
- Automatic configuration** The autosampler connected to the local area network (LAN) is automatically detected by the chromatography software.
- Device Status** When used in a local area network (LAN), the system status of the autosampler can be checked using the chromatography software.

Connecting the Device in a Local Area Network (LAN) to a Computer

In case no connection between the computer and the devices can be established via LAN, go through the following steps. Check after each step if the problem is solved. If the problem cannot be located, call the Technical Support.

<p>1. Check the status of the LAN connection in the Windows taskbar:</p> <ul style="list-style-type: none"> ▪  Connected ▪  Connection not established <p>If no connection was established, test the following:</p> <ul style="list-style-type: none"> ▪ Is the router switched on? ▪ Is the patch cable connected correctly to the router and the computer? 	<input type="checkbox"/>
<p>2. Check the router settings:</p> <ul style="list-style-type: none"> ▪ Is the router set to DHCP server? ▪ Is the IP address range sufficient for all the connected devices? 	<input type="checkbox"/>
<p>3. Check all connections:</p> <ul style="list-style-type: none"> ▪ Is the patch cable connected to the LAN ports and not the WAN port? ▪ Are all cable connections between devices and router correct? ▪ Are the cables plugged in tightly? 	<input type="checkbox"/>
<p>4. If the router is integrated into a company network, pull out the patch cable from the WAN port.</p> <ul style="list-style-type: none"> ▪ Can the devices communicate with the computer, even though the router is disconnected from the company network? 	<input type="checkbox"/>

5. In case you own a Control Unit, check the settings in the menu <i>Setup > Network</i> . <ul style="list-style-type: none">▪ Is <i>LAN-DHCP</i> set for controlling?▪ Did the device receive an IP address?	<input type="checkbox"/>
6. Turn off all devices, router, and computer. Firstly, turn on the router and secondly turn on the devices and the computer. <ul style="list-style-type: none">▪ Has this been successful?	<input type="checkbox"/>
7. Replace the patch cable to the device which no connection could be established with. <ul style="list-style-type: none">▪ Has this been successful?	<input type="checkbox"/>
8. Ensure that the IP port of the device corresponds to the IP port set in the chromatography software.	<input type="checkbox"/>

Injection Principles

The autosampler can work with the following injection principles:

ILD™ of Spark Holland

For high-pressure injections of up to 1000 bar, the autosampler is equipped with a so-called ILD™ valve (Intermediate Loop Decompression). This valve consists of a rotor-stator combination and includes a central port for depressurizing. For high-pressure applications, the sample loop is depressurized prior to receiving the sample. This way, the sample is not diluted with solvent. Because the valve is switched extremely fast, pressure spikes are reduced. Analyses are more precise and wear of the column is reduced.

Legend

- ① Sample loop
- ② Syringe
- ③ Sample vial
- ④ Column
- ⑤ Pump


Fig. 10 ILD™ injection principle

PASA™ Loop Injection Principle

The autosampler uses loop injection with pressure assistance (Pressure Assisted Sample Aspiration PASA™) and has the following features:

- Samples do not have to be degassed
- No air bubbles in sample loop
- No clogging or contamination of sample needle
- Precise control of syringe movement

Injection Methods

Three different injection methods are available:

- Full loop injection
- Partial loop injection
- Microliter pick-up injection

Full Loop Injection

In full loop injection mode, the sample loop is completely filled with the sample. The maximum reproducibility but not the maximum precision is achieved because the size of the sample

loop may have a deviation of $\pm 10\%$. The maximum injection volume equals the loop volume. The sample loop is filled with a multiple of the loop volume:

- 3 x loop volume for loops up to 100 μl
- 2 x loop volume for loops between 100 and 500 μl
- 1,5 x loop volume for loops over 500 μl

The sample loss per injection is the sum of the overfilling of the sample injection times x and the flush volume set for the needle used.

Partial loop injection

In partial loop fill mode, the sample loop is filled with both sample and solvent. This ensures the highest precision of the sample volume with minimal loss of sample. The maximum injection volume equals 50 % of the loop volume. The sample loss per injection equals the flush volume plus 3 times the sample volume set for the needle used.

Microliter pick-up injection

In microliter pick-up injection mode, the sample loop is filled with a very small amount of sample and solvent (mobile phase). This ensures very high precision with no loss of sample. The maximum sample volume is at 50 % of the loop volume minus 3 times the needle volume.



Air needle and sample needle



The autosampler uses a system of two telescopic needles, one that pierces through the cap of the sample vial – the air needle – and one that extracts the sample – the sample needle.

Using a syringe, the sample is aspirated through the two needles out of the sample vial while under pressure and into the sample loop. To prevent the syringe from becoming contaminated, a buffer tube is situated between the syringe and the valve. Using washing solution, sample residue is removed from the sample needle and buffer tube.

Full Loop Injection

Note A flushing process takes places after every injection.

Functional schematic	Explanations
 <p>Fig. 11 Full loop injection: Initial conditions</p>	<p>1. Initial conditions: The valve switches to the INJECT position. Sample needle and air needle are inserted into the vial. The air needle creates pressure, which prevents air and steam bubbles from developing.</p>
 <p>Fig. 12 Full loop injection: The needle and sample lines are flushed</p>	<p>2. The syringe aspirates the flush volume from the sample vial into the sample line and removes any washing solution.</p>

Functional schematic	Explanations
 <p>Fig. 13 Full loop injection: Valve switches to <i>LOAD</i> position</p>	<p>3. The valve switches to <i>LOAD</i> position in order to transport the sample material to the inlet of the sample loop.</p>
 <p>Fig. 14 Full loop injection: The sample loop is fully filled.</p>	<p>4. The sample loop is filled by transporting a certain amount of the loop volume (depending on the volume of the loop) through the loop. 3 x loop volume for loops up to 100 µl</p>


Functional schematic	Explanations
 <p>The diagram illustrates the full loop injection process. A syringe is connected to a buffer tubing system. A valve, shown in a circular inset with positions 1 through 6, is currently in the 'INJECT' position (position 1). This position routes the sample from the injection loop through the column. A separate line shows the needle in a 'Wash/Transport position' over a tray of vials.</p>	<p>5. The valve switches to the INJECT position. From that moment, the sample loop is part of the HPLC flow of the mobile phase. The sample is transported to the column. The analysis starts.</p>

Fig. 15 Full loop injection: Valve switches to INJECT position

Next steps Flush the needle after every injection.

Air Segment for Full Loop Injection

You can use an air segment of 5 µl to reduce the flush volume. The air segment precedes the flush segment and is not injected.

The flush volume must be at 30 µl in case of a standard needle and an injection with air segment; in case of an injection without air segment, the flush volume has to be 35 µl. You may need to insert a higher flush volume for extremely viscous samples to reduce the syringe speed and improve the performance.



Fig. 16 Air Segment for full loop injection

Partial Loop Fill

Flush volume An air segment of 5 µl preceding the flushing volume, which is not injected, can reduce the amount of sample dilution caused by the dispersion during aspiration.


- Choose the following flush volumes for a standard needle:
 - At least 30 µl for injections with air segment
 - 35 µl for injections without air segment
- Increase the flush volume and reduce the syringe speed to achieve better results for highly viscous samples.

Sample Volume In the Autosampler 3950, the syringe moves the sample into the sample loop. In partial loop fill mode, the sample volume is maximum 50 % of the loop volume.

Legend

- ① Flushing solution
- ② Tube connector for waste
- ③ Syringe
- ④ Tube
- ⑤ Valve
- ⑥ Sample loop
- ⑦ Column
- ⑧ Needle
- ⑨ Microtiter plate with samples
- ⑩ Capillary

Note The partial loop fill is processed automatically.







Figure	Process
 <p>Fig. 19 Partial loop fill: Valve position "Inject"</p>	<p>1. The valve switches to the <i>INJECT</i> position. The tube and the capillary are filled with flushing solution.</p>
 <p>Fig. 20 Partial loop fill: Aspirating the sample</p>	<p>2. The needle moves into the sample and the syringe aspirates the sample. The tube is partially and the capillary fully filled with sample.</p>
 <p>Fig. 21 Partial loop fill: Valve position "Load"</p>	<p>3. The valve switches into the <i>LOAD</i> position and the syringe aspirates the sample into the sample loop. The tube is full of flushing solution and sample, the capillary of sample and solvent.</p>

Figure	Process
 <p>Fig. 22 Partial loop fill: Injecting into the column</p>	<p>4. The valve switches to the <i>INJECT</i> position and the pump transports the sample and the solvent to the column.</p>
 <p>Fig. 23 Partial loop fill: Aspirating flushing solution</p>	<p>5. The syringe aspirates the flushing solution from the bottle.</p>
 <p>Fig. 24 Partial loop fill: Flushing sample and solvent</p>	<p>6. The syringe injects the flushing solution into tube and capillary. Both sample and solvent get flushed.</p>

Microliter Pick-Up Injection



Figure	Process
 <p>Fig. 25 Microliter pick-up injection: Initial conditions</p>	<ol style="list-style-type: none"> 1. Initial conditions: The valve switches to the INJECT position. The sample needle is in position wash/transport. The transport container is filled with washing liquid. 2. The sample line is filled with washing fluid until the fluid reaches the sample loop inlet. The valve remains in the INJECT position during filling/transportation. <p>Note: Make sure that washing liquid and solvent are compatible.</p>
 <p>Fig. 26 Microliter pick-up injection: Sample line is filled with solvent</p>	<ol style="list-style-type: none"> 3. The valve switches to LOAD. A washing fluid segment is aspirated into the sample loop.





Figure	Process
 <p data-bbox="204 840 798 873">Fig. 27 Microliter pick-up: Sample is aspirated</p>	<p data-bbox="1005 257 1396 358">4. The needle moves from the transport position to the sample vial.</p>
 <p data-bbox="204 1467 981 1500">Fig. 28 Microliter pick-up: The injection volume is aspirated.</p>	<p data-bbox="1005 896 1396 1064">5. The sample is aspirated from the sample vial according to the amount of programmed injection volume.</p>

Figure	Process
 <p data-bbox="204 817 885 884">Fig. 29 Microliter pick-up: The sample is transported through the sample loop.</p>	<p data-bbox="1011 253 1380 459">6. The needle moves back into transport position. A second washing liquid sample is aspirated. The sample is transported through the sample loop.</p>
 <p data-bbox="204 1467 965 1534">Fig. 30 Microliter pick-up: The sample is transported to the column.</p>	<p data-bbox="1011 902 1380 1142">7. The valve switches to INJECT. At that moment, the sample loop is part of the HPLC flow of the mobile phase. The sample is transported to the column. The analysis starts.</p>


Air Segment with μ l-Pick-Up Injection

If an air segment has been programmed, it will be inserted in front of the solvent area and in front of each sample.

The following conditions apply:



- The air segment is injected into the HPLC system in front of the sample segment.




- In this mode, the needle does not create pressure because the sample volume could be distorted due to the air expansion during the movement from sample vial to wash position.

Microliter Pick-Up Injection with 84+3

If you choose to operate the microliter pick-up with the 84+3 vial plate, three 10 ml vials are used automatically for transportation. In this case, the needle position wash/transport is only used for washing the needle. The following injection process is the result:

Process	Figure
1. The sample needle begins with the wash/transport position. The valve starts with the INJECT position.	 <p>Fig. 32 Microliter pick-up 84+3: Initial conditions</p>
2. The first injection starts in the transport position with the syringe aspirating transport liquid from the vial to fill the sample line with transport liquid and to remove the flushing solution.	 <p>Fig. 33 Microliter pick-up 84+3: Sample line is filled with flushing solution.</p>

Process	Figure
<p>3. The valve switches to LOAD. The programmed injection volume is aspirated out of the sample vial.</p>	 <p>Fig. 34 Microliter pick-up 84+3: Injection valve switches to <i>LOAD</i> position</p>
<p>4. The needle moves back into the transport position. A second washing liquid segment is aspirated. The sample is transported through the sample loop.</p>	 <p>Fig. 35 Microliter pick-up 84+3: The sample is transported through the sample loop.</p>
<p>5. The valve switches to INJECT. From that moment, the sample loop is part of the HPLC flow of the mobile phase. The sample is transported to the column. The analysis starts.</p>	 <p>Fig. 36 Microliter pick-up 84+3: Injection valve switches to position <i>INJECT</i>.</p>


Air Segment with μ l-Pick-Up Injection 84+3

If an air segment has been programmed, it is inserted in front of the transportation liquid segment and in front of each sample.

Details on this injection mode:

- The air segment is injected into the HPLC system in front of the sample segment

- In this mode, the needle does not create pressure because the sample volume could be distorted due to the air expansion during the movement from sample vial to wash position

84+3 Vial Plate

The 84+3 vial plate is to be ordered separately. Position the plate onto the plate holders inside the autosampler.

Hinweis: Position 87 of the vial plate must always be located at the rear right corner of the sample extraction chamber.

The vial plate offers space for maximum 84 vials with 1.5 ml and 3 vials with 10 ml. It was designed with the caps of all vials being level, independently from their unique heights, which means that one needle can be used for all vial sizes.


Fig. 38 Height of 84+3 vial plate

The length of the sample needle must be programmed at 1.5 ml. Relating to the liquid level in 10 ml vials, it is possible to program the needle at two different lengths.


Fig. 39 Needle length and liquid level

Follow the numeric order if you are to program a sample sequence for a vial plate.


Fig. 40 Sequence

In case you are creating a sequence chart with only one sample per row (1 sample in each run), programming the sequence is free to you.

µl-Pick-Up Parameters for the 84+3 Vial Plate

Before using the 84+3 vial plate, it is necessary to change the software settings.

Note If you choose to operate the microliter pick-up with the 84+3 vial plate, three 10 ml vials are used automatically for transportation. In this case, the needle position wash/transport is only used for washing the needle.

The vial positions for the 84+3 vial plate are as follows:

Position of the first sample:	Vial positions 1-84
Position of the final sample:	Vial positions 1-84
First destination position:	Vial positions 1-84
Transport position:	Vial positions 85-87


Fig. 41 Reagent/Transport position on 84+3 vial plate

You need to define the first and the final sample to enter a sample range. The transport positions are always positions 85, 86, 87.

Details on Programming the 84+3 Injection Method

- The position of the transport vial can be programmed. Some possible positions are 85, 86, and 87. Position 85 is the standard position for transport vials. The system calculates the necessary transport volume. Make sure before system start that the vial is filled with at least 8000 µl of liquid.
- The liquid levels of the transportation liquids are not updated at system start. To avoid contamination of the air needle, the needle stops inside the transport vial in the up most position.

- If conducting a sample sequence or one sample per row, the autosampler takes the residual volume of the transport liquid into account. If the volume falls below 4000 µl, the needle moves deeper into the transport vial. If the amount of transport liquid reaches 0 µl, the autosampler sends out error message 369 (not enough transport liquid available).
- The liquid levels of the transportation liquid will reset after reprogramming the mode.
- The needle does not move automatically to the next transport vial. If transport liquid is to be taken from another position, you have to change the program settings.

Air Needles

Six different lengths of air needles from 50 to 80 mm are available for the autosampler. The needle holder allows you to further adjust the needle height by 6 mm.

Standard Air Needle

The standard air needle is 62 mm long and can be used for a wide range of deep and shallow vial plates.

When 10 ml sample vials are used, the needle deeply penetrates the sample vial. If this is not filled to more than 60%, the needle can be used in the typical manner. The same applies to deep microtiter plates.

For non-standard settings, use the corresponding needle types.

Legend

- ① 10 ml sample vial
- ② 2 ml sample vial


Fig. 42 Standard air needle with 10 and 2 ml sample vial

Note The PASA™ loop injection principle is not suitable for shallow microtiter plates. The function of the air needle is only ensured when it pierces the closure to a sufficient degree.

Legend

- ① Deep microtiter plate with closure.
- ② Shallow microtiter plate


Fig. 43 Standard air needle with microtiter plates

Choosing the Correct Air Needle

To choose the correct air needle, take the following dimensions into consideration:

Legend

H_t = height of vial plate

D_w = hole depth

C_d = thickness of closure

N_h = set needle height

A_c = distance from air needle tip to closure (min. 2 mm)

? = excess length


Fig. 44 Calculating the correct air needle

Condition: $H_t - D_w = 2$ to 6 mm

Excess length of the sample needle:

$$H_t - C_d - N_h - A_c = ?$$

Choose the correct needle type on the basis of the excess length.

Air needle type	Protrusion length
50 mm, yellow	34-40 mm
56 mm, red	28-34 mm
62 mm, white (standard needle)	22-28 mm
68 mm, blue	16-22 mm
74 mm, green	10-16 mm
80 mm, black	4-10 mm

Legend

- ① 10 ml sample vial,
50 mm air needle
- ② 2 ml sample vial,
62 mm air needle


Fig. 45 Selected air needles with sample vials

Legend

- ① Deep microtiter plate with closure.,
56 mm air needle
- ② Shallow microtiter plate,
80 mm air needle


Fig. 46 Selected air needles with microtiter plates

Calculation Example for Air Needle

Starting point:

- Autosampler with standard setting for needle height.
- Deep microtiter plate with closure.

Dimensions:

$$H_t = 41.4 \text{ mm}$$

$$D_w = 37.8 \text{ mm}$$

$$C_d = 3.8 \text{ mm}$$

$$N_h = 6.0 \text{ mm (standard)}$$

$$A_c = 2.0 \text{ mm (minimum)}$$

Condition:

$$H_t - D_w = \text{between 2 and 6 mm}$$

$H_t - D_w = 41.4 \text{ mm} - 37.8 \text{ mm} = 3.6 \text{ mm}$. Condition has been met.

$$H_t - C_d - N_h - A_c = \text{excess length}$$

$$41.4 \text{ mm} - 3.8 \text{ mm} - 6.0 \text{ mm} - 2.0 \text{ mm} = 29.6 \text{ mm}$$

Air needle type	Protrusion length
56 mm, red	28-34 mm

An air needle length of 56 mm is required.

Handling the Sample Vials

When handling the sample vials, observe the following:

- Fill the sample vials using a pipette to allow air to escape.
- To prevent the sample from contaminating the air needle, do not fill the sample vials to the very top.
- Do not use sample vials that are unclosed.
- Only use air-tight closure seals to prevent air bubbles from forming and volatile components from evaporating.
- Do not use sample vials with hard closures that the sample needle cannot pierce.

Mixing and Diluting

A mix method can be programmed for the autosampler to mix or dilute the sample fluid.

- Configure the mixing routine and syringe speeds using the chromatography software.
- A maximum of 15 steps can be programmed for a mix method.

Three types of actions are possible:

1. (*Add*)
2. (*Mix*)
3. (*Wait*)

Add When adding, the defined volume is aspirated from either the sample vial, the vial with *Reagent A* or *Reagent B* or flushing fluid and then dispensed into the destination vial.

Note To prevent carryover, the autosampler removes 125% of the given volume from the corresponding sample vial and uses the additional 25% to flush the tube and needle.

Mix When mixing, the contents of a specific sample vial is mixed by aspirating and dispensing the defined volume *n* times. If a destination vial has not been defined, mixing is performed in the current sample vial.

Note When entering the *Sample Vials*, the *Destination Vial* is used automatically.

Wait With the *Wait* command, the system waits until the programmed delay time has elapsed before executing the next line of the program.

An example: Add

The command *ADD 100 µl from Reagent A to Destination* involves the following steps:

1. An air segment of 5 µl is aspirated to separate the flushing solution in the buffer tube from *Reagent A*.
2. 50 µl of *Reagent A* are aspirated to flush the tube and needle.
3. Syringe is emptied into the waste container through the drainage tube.
4. 100 µl of *Reagent A* are aspirated and then dispensed into the destination vial.
5. Tube and needle are flushed with flushing solution.

An example: Mix

If this is preceded by an *ADD... to Destination* command, mixing is performed in the destination vial. If this is preceded by an *ADD to Sample* command, mixing is performed in the sample vial.

The *MIX 3 times with 100 µl* command triggers the following steps:

1. An air segment of 50 µl is aspirated to separate the flushing solution in the buffer tube from the sample solution to be mixed.
2. Syringe is emptied into the waste container through the drainage tube.
3. 100 µl solution are aspirated and dispensed into the same sample vial.
4. Step 3 is repeated twice.
5. Tube and needle are flushed with flushing solution.

Sample Positions in Mix Method

When configuring a mix method, the positions of the sample vials depend on whether the vial plates are to be processed in rows or columns.

Processing in Columns

When column processing is used, the following positions are possible for the *Sample*, *Destination*, *Reagent A* and *Reagent B*:



Fig. 47 Sample positions in columns

Processing in Rows

When row processing is used, the following positions are possible for the *Sample*, *Destination*, *Reagent A* and *Reagent B*:

Parameters for Mixing Method with 84+3 Vial Plate

Before using the 84+3 vial plate, it is necessary to change the software settings.

The vial positions for the 84+3 vial plate are as follows:

Position of the first sample:	Vial positions 1-84
Position of the final sample:	Vial positions 1-84
First destination position:	Vial positions 1-84
Reagent position:	Vial positions 85-87


Fig. 49 Reagent/Transport position on 84+3 vial plate

You need to define the first and the final sample to enter a sample range. The same range can be used for the destination vials. The reagent positions are always positions 85, 86, 87.

Details on Programming the 84+3 Mixing Method

- The position of the reagent vial can be programmed. Some possible positions are 85, 86, and 87. The standard transport positions are 86 for Reagent A and 87 for Reagent B. The system calculates the necessary reagent volume. Make sure before system start that the vials are filled with at least 8000 μl of liquid.
- The liquid levels of the reagent liquids are not updated at system start. To avoid contamination of the air needle, the needle stops inside the reagent vial in the up most position.
- If conducting a sample sequence or one sample per row, the autosampler takes the residual volume of the reagent liquid into account. If the volume falls below 4000 μl , the needle moves deeper into the reagent vial. If the amount of reagent liquid reaches 0 μl , the autosampler sends out error message 370 (not enough reagent liquid available).
- The liquid levels of the reagent liquid will reset after reprogramming the mode.
- Reagent will only be taken from preprogrammed positions. The needle does not move automatically to the next reagent vial. You have to change the settings when reagent is to be taken from one of the other vials.

Capillary and Hose Connections

Note Flush the system before connecting the column.

The connection scheme inside the autosampler displays the correct connections:

Legend

- ① Waste hose
- ② Flushing fluid container
- ③ Buffer tube
- ④ Sample loop
- ⑤ Column
- ⑥ Pump
- ⑦ Injection needle
- ⑧ Position: Injection
- ⑨ Position: Load


Fig. 50 Connection scheme for capillaries and tubing

Connecting the Injection Valve

Legend

- ① Connection tube buffer to syringe
- ② Drainage for ILDTM valve
- ③ Connection plastic capillary to injection needle
- ④ Connection sample loop
- ⑤ Connection sample loop
- ⑥ Connection stainless steel capillary to column
- ⑦ Connection stainless steel capillary to pump


Fig. 51 Connections on the injection valve

Connecting the Syringe

Legend

- ① Connection transport liquid
- ② Connection buffer tube
- ③ Connection flushing solution



Fig. 52 Syringe connections

Tube Guide for Flushing Solution

Note Use the tube guide in the collecting container for the flushing solution in order to not hinder horizontal movement of the needle.

Legend

- ① Tube Guide for Flushing Solution
- ② Collecting container with hole for tube guide (flushing fluid)

**Fig. 53 Tube guide for flushing solution tube**

Connecting the Drainage Tubing

The waste drainage removes all flushing fluids and non-injected sample solutions.

- Connect the condensed water and drainage tubing to the front panel of the device.
- Connect the adapter and insert the end of the tube into a container on the floor.

Note Make sure that the drainage tubing is not pinched so that the liquid can drain.

Legend

- ① Tube for waste liquid
- ② Tube for condensed water
- ③ Adapter for drainage tube


Fig. 54 Connecting the drainage tubing

PEEK Connections

Disposable PEEK fittings

- One-piece disposable polyetheretherketone fittings (PEEK) for easier installation of flexible, thin capillaries (outer diameter 0.5 mm)
- Tightening torque of the PEEK screw:
Tightened by hand (about 0.5 Nm)

Connecting the Autosampler to Other Devices

Controlling the Autosampler with Chromatography Software

You control the autosampler directly with software, e. g. KNAUER OpenLAB[®] and ClarityChrom[®].

Using the LAN connection at the rear side of the device, connect the autosampler to the network.

Checking and Configuring the Parameters of the Autosampler

You can configure the autosampler parameters using the chromatography software, e.g. ClarityChrom[®]:

1. Select the autosampler in LAN.
2. Set the syringe volume to either 250 μl (default) or 500 μl .
3. Enter the serial number of the autosampler.
4. Choose the cooling option when your autosampler offers temperature control.

Configuration Window of ClarityChrom[®]

Legend

- ① Serial number
- ② Device detection in local network
- ③ Manual search for device in network
- ④ Volume of sample loop
- ⑤ Volume of syringe


Fig. 55 Autosampler configuration with ClarityChrom[®]

Autosampler Device Software

A software CD with *Autosampler 3950 Service Manager* is included with the autosampler. Using this software, you can check and control the device. These settings can also be made with the chromatography software, so that it is not absolutely necessary to install the *Service Manager*.

Flushing the System

The system should be flushed before the column is connected. The flushing of the system can be controlled using the chromatography software or the *Autosampler 3950 Service Manager*.

Note KNAUER recommends using a mixture of water and isopropanol (80%/20%) or the mobile phase as the flushing solution.

The following steps are explained for when the *Service Manager* is used.

System Flushing with *Autosampler 3950 Service Manager*

1. Install *Autosampler 3950 Service Manager*.
2. Fill the flushing solution into a solvent bottle and degas it using helium or an ultrasonic bath.
3. Slide the tube for the flushing solution into the solvent bottle.
4. Select the *Alias* \triangleright \Rightarrow *Direct Control* menu.
5. In the *Syringe* field, click *End*. One syringe volume is aspirated into the syringe through the flushing solution tube.
6. In the *Syringe* field, click *Home*. The syringe content is emptied into the drainage tube.
7. Repeat step 5 and 6 until the syringe and the flushing solution tube have been filled completely.
8. In the *Initial wash* field, click *Start*. All tubes that are connected to the syringe are flushed.
9. In the *Initial wash* field, click *Stop*.
10. Click *Close* to exit the *Direct Control* window.


Fig. 56 System flushing with the *Service Manager* software

I/O Connection

By default, the autosampler has an I/O connection that supports TTL inputs (low-active) and a closed-contact output. Devices without LAN connection that require a trigger signal for the injection can be connected by means of the I/O connection.

The TTL inputs allow you to control the autosampler using other devices and are defined using the instrument methods of the chromatography software.

Note Only connect the autosampler to devices that fulfill the required safety standards!

Defining the TTL Inputs

- *Next Injection Input*: The injection sequence is started. After ending the injection sequence, the autosampler waits for the next start signal.
- *Freeze Input*: The analysis time is paused. The autosampler runs the configured program up to the filling of the sample loop. The injection is not performed until the input is deactivated.
- *Stop Input*: Immediately halts the analysis.

Defining the Closed-Contact Output

- *Inject Marker*: The closed-contact output is activated when the injection valve switches from *LOAD* to *INJECT*.
- *Alarm*: The closed-contact output is activated when there is an autosampler fault.
- Auxiliary: -

Configuration of I/O Connection (9 pins)

Explanation	Cable color
1. Output, start injection	Red in three-conductor cable
2. Output, start injection	Black in three-conductor cable
3. Input 1, programmable input for stopping injections (low-active)	Red in four-conductor cable
4. Input 2, programmable input for stopping injections (low-active)	Black in four-conductor cable
5. Not occupied.	-
6. Output	Brown in three-conductor cable.
7. Not occupied.	-
8. Ground, for inputs 1 and 2	Orange in four-conductor cable
9. Ground, for inputs 1 and 2	Brown in four-conductor cable

Device Test

The reproducibility of the sample volume is a critical factor for maintaining high-quality analysis results.

- Test the autosampler features with the chromatography software on a regular basis.
- If the device test determines that the autosampler does not fulfill the requirements, mark the device as defective and do not continue using it.
- Do not re-use the autosampler until it has been repaired and/or serviced.

Test Intervals

Run the device test at the following time intervals:

- Average use of 1 to 5 days/week: Device test every 6 months
- Average use of more than 5 days/week or 24 hours/day: Device test every 3 months
- Operation with buffer solutions or other salt solutions: Device test every 3 months

Devices and Components for the Test

- Autosampler 3950, standard version with 10 µl sample loop and 250 µl syringe
- HPLC pump, 1 ml/min flow rate
- UV detector (with a data rate of 50 Hz if possible, otherwise 10 Hz, flow cell: 10 mm path length)
- Chromatography software
- Solvent: 90 % water, 10 % methanol (HPLC quality)
- Test solution (sample):
 - *a*: 50 ppm uracil dissolved in water (HPLC quality)
 - *b*: 250 ppm uracil dissolved in water (HPLC quality)
- Flushing solution:
 - 80 % water, 20 % isopropanol (HPLC quality)
 - Alternative: 80 % water, 20 % methanol (HPLC quality)
- Restriction capillary:
 - Inner diameter 0.25 mm
 - Length 200 cm

Note Degas the solvent to prevent malfunctions caused by the presence of air bubbles.

1. Reproducibility of Sample Volume

The variation coefficient must not exceed 0.5 %.

- Configure the HPLC pump, UV detector and autosampler using the chromatography software
- Inject 10 µl of test solution *a* (50 ppm uracil, dissolved in water)

Standard Setting of Autosampler

- Loop volume: 10 µl
- Tubing volume: 15 µl
- Syringe volume: 250 µl

Method Parameters of Pump

- Flow: 1 ml/min
- Time: 1 min

Method Parameters of Autosampler

- Injection method: Partial loopfill
- Syringe speed: normal
- Flush volume: 30 µl, (40 µl for microtiter plates)
- Needle wash: active 2 times
- Air segment: yes
- Headspace pressure: yes
- Injections/vial: 9 (7 for microtiter plates)
- Inj. volume: 2 µl
- Vial position: 1A1. Start the single run with

Method Parameters of UV Detector

- Wavelength: 254 nm
- Sampling rate: 50 Hz if possible, else 10 Hz
- Time: 0.5 min

Configuring Repeat Runs of Autosampler

- Injections/vial: 9x (7x for microtiter plates)
- Inj. volume: 1 µl

Starting Repeat Runs

1. Put a vial with at least 500 µl test solution on position *1A1* of the sample plate.
2. Start the repeat runs.

Analyzing the Individual Chromatograms

1. Calculate the average of the measuring values of the *Peak areas*.
2. Calculate the variation coefficient VK_1 .
3. Enter the results into the *Test Report* form.

Formula for Determining the Arithmetic Mean

$$\overline{Peakareas} = \frac{\sum_{i=1}^n Peakareas_i}{n}$$

- Formula for determining the standard deviation ($i = 1-9$):

$$\sigma_{n-1} = \sqrt{\frac{\sum (Peakarea_i - \overline{Peakareas})^2}{n-1}}$$

- Formula for determining the variation coefficient:

$$VK_1 = \frac{\sigma_{n-1}}{\overline{peakarea}} \times 100\%$$

2. Sample Carryover

The percentage of sample carryover must not exceed 0.3 %.

- Alternately inject 10 μ l of test solution *b* (250 ppm uracil, dissolved in water) followed by solvent.
- Position of sample vial: *1A1*
- Position of solvent: *1A2*
- Injection volume: 1 μ l

Creating a Sequence with 6 Lines

- Test solution: Position *1A1*
- Solvent: Position *1A2*
- Injection volume: 1 μ l
- Repeats: 1

Note For microtiter plates, select 6 consecutive positions that are alternately to be filled with test solution and solvent.

Analyzing the Individual Chromatograms

1. Calculate the average of the measuring values of the *Peak* areas.
2. Put the average of the solvent injection in relation to the average of the test solution injection.
3. Enter the results into the *Test Report* form.

Formula for Calculating Sample Carryover

$$PV = \frac{\sum_i \frac{Peakarea_{i\text{ Fließmittel}}}{3}}{\sum_i \frac{Peakarea_{i\text{ Testlösung}}}{3}} \times 100 \%$$

3. Linearity

To determine the linearity, the correlation coefficient of the regression lines is determined from the measured values for the *peak* areas and injection volume.

The correlation coefficient must not exceed 0.998 %.

- Inject 10, 20, 30, 40 and 50 µl of test solution *b* (250 ppm uracil, dissolved in water) respectively.
- Position of sample vial: 1A1
- Injection volume: 10 µl, 20 µl, 30 µl, 40 µl, 50 µl
- Repetition: 3

Note Fill consecutive positions on the microtiter plates with test solution.

Analyzing the Individual Chromatograms

1. Calculate the correlation coefficient *r* of the regression lines from the measured values for the *Peak* areas and the injection volume.
2. Enter the results into the *Test Report* form.

Formula for Determining the Correlation

$$r = \frac{\sum x_i y_i - n \bar{x} \bar{y}}{\sqrt{(\sum x_i^2 - n \bar{x}^2)(\sum y_i^2 - n \bar{y}^2)}}$$

y_i = Y value of measured value i (injection volume)

x_i = X value of measured value i (peak area)

\bar{y} = arithmetic mean of Y across all n measured values

\bar{x} = arithmetic mean of X across all n measured values

n = number of measured value pairs

4. Mixture Test

To create a mixing method, follow the instructions in the manual of the chromatography software.

- The test solution variation coefficient $VK2$ must not exceed 0.5 %.
- The dilution variation coefficient $VK2$ must not exceed \leq 0.5 %.
- The dilution factor $F10$ has to be within the range of $9.85 < x < 10.25$.
- Test solution: 10 μ l (50 ppm uracil, dissolved in water)
- For injecting the dilution, go to *Mix methods* and create a mixing method in which 40 μ l of the test solution is mixed with 360 μ l solvent.
 - Dilution: 10 μ l, 5 ppm uracil, dissolved in deionized water
- Inject 10 μ l test solution and und 10 μ l dilution three times each.

Creating a Sequence with 2 Lines

- Injection volume: 1 μ l
- Repeats: 3

Positioning Vials for Dilution

- To position the vials on the sample plates, select the *Columns* option in the chromatography software.
- Put a vial with test solution (*Sample*) at position *2A1* of the vial plate.
- Put an empty vial (*Destination*) at position *2A5* of the vial plate.
- Put a vial with solvent (*Reagent A*) at position *1A1* of the vial plate.

Note Observe the plate assignments for the vials (*Sample*, *Reagent A*, *Reagent B*, *Destination*).

Analyzing the Individual Chromatograms

1. Calculate the mean of the measuring values of the *Peak* areas for the test solution and the dilution.
2. Calculate the variation coefficient VK_2 .
3. Calculate the variation coefficient VK_3 .
4. Calculate the dilution factor *F10* from the ratio of the mean of the *Peak* areas of the test solution and the dilution.
5. Enter the results into the Test Report form.

Archiving

- Enter all test results into the *Test Report* form.
- Enter the serial number, date of the test, date of the next test and name of the tester.
- File the *Test Report* form in the device logbook.

Test Report

Module	Autosampler	
Autosampler 3950	Standard <input type="checkbox"/>	Sample cooling <input type="checkbox"/>
Serial number		

No	Test	Setting	Specification	Result
1	Reproducibility	<ul style="list-style-type: none"> ▪ Inject 10 µl test solution nine times. ▪ Microtiter plate: Inject 10 µl test solution seven times. 	$VK_1 \leq 0.5 \%$	
2	Carryover	Alternately inject 10 µl test solution and 10 µl solvent three times.	$PV \leq 0,3 \%$	
3	Linearity	Inject 10, 20, 30, 40 and 50 µl test solution three times each.	$r \geq 0,998$	
4	Mixture test	Inject 10 µl test solution and 10 µl of the dilution created by the autosampler, three times each.	$VK_2 \leq 0.5 \%$ $VK_3 \leq 0.5 \%$ F10: $9,85 < x < 10,25$	

Date:

Date of the next device test:

Tester:

Operation Qualification Report

Legend

- ① Detailed information of the autosampler
- ② Detailed function test
- ③ Date and fields for signature and other information


Fig. 57 Operation Qualification Report

Maintenance and Care

Contacting the Technical Support

You have various options to contact the Technical Support:

Phone +49 30 809727-111

Fax +49 30 8015010

E-mail support@knauer.net

You can make your requests in English and German.

Maintenance Contract

The following maintenance work on the device may only be performed by KNAUER or a company authorized by KNAUER and is covered by a separate maintenance contract:

- Opening the device
- Removing the hood or the side panels.

Which Type of Maintenance Tasks May Users Perform on the Device?

Users may perform the following maintenance tasks themselves:

- Exchanging the fuses
- Exchanging the air and sample needle
- Exchanging the injection valve
- Exchanging the sample loop
- Exchanging the rotor seals
- Exchanging the capillary and tubing

CAUTION!

Performing maintenance tasks on a switched on device can cause damage to the device.

Switch off the device and pull the power plug.

Leaks in the Capillary Screw Fittings

Note If leaks occur on the capillary screw fittings after maintenance and proper assembly, do not tighten them further, but instead replace them with new connection capillaries.

Exchanging the Fuses

1. Switch off the autosampler and remove the power plug to completely disconnect the device from the power supply.
2. Remove the fuses from the fuse box at the rear of the device.
3. Insert new fuses (2 x 2.5 A).

4. Plug in the power plug.




Changing the Valve and Rotor Seal

- Regularly clean the rotor seal of the injection valve.
- Regularly replace the rotor seal (approx. every three years).
- Remove the front panel of the autosampler.
- Remove the capillary connections, except the sample loop, from the valve.
- During removal, consecutively loosen all screws by half a turn respectively, until they can be removed.

Note Do not remove the screw in the cover plate hole diagonal to the valve!

Removing the Injection Valve and Rotor Seal

Procedure

Process	Figure
<ol style="list-style-type: none"> 1. With a screwdriver, remove screw ① on both sides of the injection valve housing. 2. Remove the injection valve. 3. With a hexagon wrench key, remove the screws ② from the valve stator block ③. 4. Carefully remove the stator. 5. Remove the rotor seal ④ from the rotor ⑤. 6. Clean or exchange the rotor seal. 	 <p data-bbox="922 689 1342 723">Fig. 58 Disassembling the valve</p>  <p data-bbox="922 1099 1342 1133">Fig. 59 Removing the rotor seal</p>  <p data-bbox="922 1420 1294 1453">Fig. 60 Rotor and rotor seal</p>

Installing the Injection Valve

- During installation, hold the injection valve in such a way that the bore hole for connecting the steel capillary to the pump (port 1) is facing upward.
- Alternately tighten all screws by half a turn, until all screws have been fully tightened.

Legend

- ① Labeling, port 1
- ② Connection from capillary to pump


Fig. 61 Installing the valve

- Procedure**
1. Insert the rotor seal.
 2. Place the valve stator block on the valve rotor and use a hexagon wrench key to tighten the screws.
 3. Insert the injection valve and use a screwdriver to tighten the screw on both sides of the valve housing.
 4. Connect the capillary again.

System Flushing

- Procedure**
1. Connect the autosampler to the power supply.
 2. Establish a connection to the PC.
 3. Start *Autosampler 3950 Service Manager*.
 4. Select the *Alias⇒Direct Control* menu.
 5. Click *Initialize* to check whether the valve is correctly positioned at the Inject position.
 6. In the *Initial Wash* field, click *Start* to flush the system.
 7. In the *Initial Wash* field, click *Stop* to stop flushing the system.

Exchanging the Sample Loop

By standard, the autosampler is equipped with a 10 µl sample loop.

- When assembling a sample loop with a different injection volume, make sure to use the correct combination of syringe and capillaries and configure the controller software appropriately.
- Always connect the sample loop to ports 2 and 5 of the injection valve.
- Calculate the maximum injection volume according to the following formula:
 - Full loop injection:
Maximum injection volume = 3 x loop volume for loops with 100 µl, 2 x loop volume for loops with 100–500 µl, 1.5 x loop volume for loops with over 500 µl.
 - Partial loop injection:
Maximum injection volume = 50 % loop volume
 - *Microliter pick-up* injection:
Maximum injection volume = 50 % loop volume - 3 x needle volume

Exchanging the Sample Needle

- When using sample plates with 12 or 48 sample vials, make sure that the needle height setting is >2 mm to prevent the needle from contacting the bottom of the sample vial.
- Tighten the screw fitting until it is finger-tight to prevent the plastic capillary from becoming blocked.

Legend

- ① Fitting
- ② Plastic capillary
- ③ Union nut
- ④ Sample needle


Fig. 62 Exchanging the sample needle

- Procedure**
1. Start *Autosampler 3950 Service Manager*.
 2. Select the *Alias*⇒*Direct Control* menu.
 3. In the *Needle* field, click *Exchange*. The needle moves to the replacement position.
 4. Loosen the union nut ③.
 5. Loosen the screw fitting ① of the plastic capillary ② on the injection valve.
 6. Remove the sample needle ④ with the plastic capillary.
 7. Install a new sample needle unit. Make sure that the air seal fully surrounds the sample needle.
 8. Fasten the sample needle with the union nut.
 9. Fasten the plastic capillary using the screw fitting on the injection valve.
 10. In the *Direct Control* window, click *Initialize*. The needle moves to the initial position.
 11. In the *Initial Wash* field, click *Start* to flush the system.
 12. In the *Initial Wash* field, click *Stop* to stop flushing the system.

13. Click *Close* to exit the *Direct Control* window.
14. Select the *Alias*⇒*Adjustments* menu.
15. On the *Needle-Tray* tab, update the settings for the sample plates.

Exchanging the Air Needle

- When exchanging the air needle, make sure that the thread of the new height adjustment screw is flush with the lower edge of the retaining nut.
- Make sure that the sealing ring is located in the retaining nut.

Legend

- ① Union nut
- ② Retaining nut
- ③ Height adjustment screw
- ④ Air needle
- ⑤ Sample needle


Fig. 63 Exchanging the air needle

- Procedure**
1. Start *Autosampler 3950 Service Manager*.
 2. Select the *Alias*⇒*Direct Control* menu.
 3. In the *Needle* field, click *Exchange*. The needle moves to the replacement position.
 4. Loosen the union nut ①.
 5. Loosen the screw fitting of the plastic capillary on the injection valve.
 6. Remove the sample needle ⑤ with the plastic capillary.
 7. Loosen the retaining nut ② and pull it downwards together with the air needle ④.
 8. Unscrew the retaining nut from the height adjustment screw ③.
 9. Screw a new air needle with a new height adjustment screw into the retaining nut.
 10. Screw in the retaining nut.
 11. Insert the sample needle and fasten with the union nut.


12. Fasten the plastic capillary using the screw fitting on the injection valve.
13. In the *Direct Control* window, click *Initialize*. The needle moves to the initial position.
14. In the *Initial Wash* field, click *Start* to flush the system.
15. In the *Initial Wash* field, click *Stop* to stop flushing the system.
16. Click *Close* to exit the *Direct Control* window.
17. Select the *Alias*⇒*Adjustments* menu.
18. On the *Needle-Tray* tab, update the settings for the sample plates.

Exchanging the Syringe

By standard, the autosampler is equipped with a 250 µl syringe. Use isopropanol as flushing solution to remove air bubbles from the new syringe.

Legend

- ① Syringe valve
- ② Syringe
- ③ Syringe drive
- ④ Syringe plunger

Exchanging the Syringe

1. Start *Autosampler 3950 Service Manager*.
2. Select the *Alias*⇒*Direct Control* menu.
3. In the *Syringe* field, click *Exchange*. The syringe plunger is lowered.
4. Unscrew the syringe ② by rotating it counterclockwise; leave the adapter in the syringe valve ①.
5. Remove the syringe plunger ④ from the syringe drive ③.
6. Fill new syringe with flushing solution.
7. Insert the syringe plunger into the syringe drive.
8. Tighten the syringe in the syringe valve by rotating it clockwise.

9. In the *Syringe* field, click *Home*. The syringe content is emptied into the drainage tube.
10. If there is still air in the syringe, click *End* in the *Syringe* field. One syringe volume is aspirated into the syringe through the flushing solution tube.
11. In the *Syringe* field, click *Home*. The syringe content is emptied into the drainage tube.
12. Slightly tap the body of the syringe if it still contains air. If necessary, repeat step 10 and 11. Repeat the step.
13. In the *Initial Wash* field, click *Start* to flush the system.
14. In the *Initial Wash* field, click *Stop* to stop flushing the system.
15. Click *Close* to exit the *Direct Control* window.

Exchanging the Syringe Plunger or Plunger Tip

1. Start *Autosampler 3950 Service Manager*.
2. Select the *Alias⇒Direct Control* menu.
3. In the *Syringe* field, click *Exchange*. The syringe plunger is lowered.
4. Remove the syringe (see above).
5. Pull the syringe plunger out of the glass cylinder of the syringe.
6. Use a pair of tweezers to remove the plunger tip.
7. Wet the new plunger tip with isopropanol.
8. Mount the new plunger tip onto the syringe plunger.
9. Push the syringe plunger into the glass cylinder of the syringe.
10. Install the syringe (see above.)
11. In the *Syringe* field, click *Home*. The syringe content is emptied into the drainage tube.

Exchanging the Syringe Valve

The syringe valve has four connections, one of them remains unused.

- Hand-tighten all fittings on the connections to the syringe valve.
- To exchange the valve, set it to the Waste position, because the hexagon socket screws are only accessible if the valve is in this position.

Legend

- ① Upper hexagon socket screw
- ② Lower hexagon socket screw
- ③ Connection flushing solution tube (hidden)
- ④ Connection buffer tube
- ⑤ Connection syringe
- ⑥ Unused connection


Fig. 65 Exchanging the syringe valve

- Procedure**
1. Start *Autosampler 3950 Service Manager*.
 2. Select the *Alias⇒Direct Control* menu.
 3. In the *Syringe* field, click *Exchange*. The syringe plunger is lowered.
 4. Loosen the lower hexagon socket screw by 2 turns.
 5. Loosen the upper hexagon socket screw by 2 turns.
 6. Pull out the upper part of the syringe.
 7. Remove the syringe.
 8. Exchange the syringe valve.
 9. Insert a new syringe.
 10. Tighten the hexagon socket screws.

Cleaning and Caring for the Device

CAUTION!

**Intruding liquids can cause damage to the device!
Place solvent bottles next to the device or in a solvent tray.**

Moisten the cleaning cloth only slightly.

All smooth surfaces of the device can be cleaned with a mild, commercially available cleaning solution, or with isopropanol.

- Clean collecting container and vial plates with a soft cloth.
- To remove deposits, flush the drainage tube regularly using solvent.

Disposal

	Hand in old devices or disassembled old components at a certified waste facility, where they will be disposed of properly.
AVV marking in Germany	According to the German "Abfallverzeichnisverordnung" (AVV) (January, 2001), old devices manufactured by KNAUER are marked as waste electrical and electronic equipment: 160214
WEEE registration	KNAUER as a company is registered by the WEEE number DE 34642789 in the German "Elektroaltgeräteregister" (EAR). The number belongs to category 8. All distributors and importers are responsible for the disposal of old devices, as defined by the WEEE directive . End-users can send their old devices manufactured by KNAUER back to the distributor, the importer, or the company free of charge, but would be charged for the disposal.
Solvents and other operating materials	All solvents and other operating materials must be collected separately and disposed of properly. All wetted components of a device, e. g. flow cells of detectors or pump heads and pressure sensors for pumps, have to be flushed with isopropanol first and water afterwards before being maintained, disassembled or disposed.

Troubleshooting

Device Errors

One possible cause of device errors is a malfunctioning valve.

Checking the Valve

Remove the valve and check all parts for wear and contamination. After the problem has been eliminated and the valve reinstalled, perform the following steps:

- Procedure**
1. Select the *Alias⇒Direct Control* menu.
 2. In the *Direct Control* window, click *Initialize*. The needle moves to the initial position.
 3. In the *Initial Wash* field, click *Start* to flush the system.
 4. In the *Initial Wash* field, click *Stop* to stop flushing the system.
 5. Click *Close* to exit the *Direct Control* window.

LAN

Software faults can occur due to flawed communications between the devices or incorrect installation of the software.

- Procedure**
1. Check the cable connections.
 2. Start *Autosampler 3950 Service Manager*.
 3. Select the *Alias⇒Direct Control* menu.
 4. In the *Direct Control* window, click *Initialize*.

In case no LAN connection between the computer and the devices can be established, go through the following steps. Check after each step if the problem is solved. If the problem cannot be located, call the Technical Support.

<p>1. Check the status of the LAN connection in the Windows task bar:</p> <ul style="list-style-type: none"> ▪  Connected ▪  Connection not established <p>If no connection was established, test the following:</p> <ul style="list-style-type: none"> ▪ Is the router switched on? ▪ Is the patch cable connected correctly to the router and the computer? 	<input type="checkbox"/>
<p>2. Check the router settings:</p> <ul style="list-style-type: none"> ▪ Is the router set to DHCP server? ▪ Is the IP address range sufficient for all the connected devices? 	<input type="checkbox"/>

<p>3. Check all connections:</p> <ul style="list-style-type: none"> ▪ Are the patch cables connected to the LAN ports and not the WAN port? ▪ Are all cable connections between devices and router correct? ▪ Are the cables plugged in tightly? 	<input type="checkbox"/>
<p>4. If the router is integrated into a company network, pull out the patch cable from the WAN port.</p> <ul style="list-style-type: none"> ▪ Can the devices communicate with the computer, even though the router is disconnected from the company network? 	<input type="checkbox"/>
<p>5. In case you own a Control Unit, check the settings in the menu <i>Setup > Network</i>.</p> <ul style="list-style-type: none"> ▪ Is <i>LAN-DHCP</i> set for controlling? ▪ Did the device receive an IP address? 	<input type="checkbox"/>
<p>6. Turn off all devices, router, and computer. Firstly, turn on the router and secondly turn on the devices and the computer.</p> <ul style="list-style-type: none"> ▪ Has this been successful? 	<input type="checkbox"/>
<p>7. Replace the patch cable to the device that no connection could be established with.</p> <ul style="list-style-type: none"> ▪ Has this been successful? 	<input type="checkbox"/>
<p>8. Ensure that the IP port of the device corresponds to the IP port set in the chromatography software.</p>	<input type="checkbox"/>

Analytical Errors

Possible causes:

- Wear due to errors in the injection and method settings.
- Unsuitable combination of sample loop, buffer tube and syringe.
- External effects such as temperature, and light-sensitive samples being exposed to light.

Solutions:

- Check whether the application has run previously without errors and that no changes have been made to the analytical system.
- Determine whether the fault is caused by the autosampler or other devices in the system.

If the required degree of reproducibility is not achieved, check the following possible sources of error and instigate steps to eliminate them:

Cause of fault	Elimination
Air in liquid path	Initialize the Smartline Autosampler 3950.
Leaking syringe	<ul style="list-style-type: none"> ▪ If the syringe is leaking at the top, check whether it has been installed correctly. ▪ If the syringe is leaking at the bottom, exchange the syringe plunger.
Leaking syringe valve	Check valve and exchange if required.
Rotor seal worn	Exchange the rotor seal and check the stator block of the valve.
Dead volume in capillary connections	Install new fittings onto capillary connections.

If an empty sample run returns an excessively large peak, check the following possible causes of error and instigate steps to eliminate them:

Cause of fault	Elimination
Solubility problems	Either modify sample or accept carry-over.
Interaction between the empty sample and the <i>hardware</i>	<ul style="list-style-type: none"> ▪ Check <i>hardware</i>: <ul style="list-style-type: none"> ▪ Flush needle (inside and outside) or install a different needle type (steel, PEEK or with glass coating). ▪ Valve: Exchange rotor seal (other material). ▪ Capillaries and tubing: Use other connections between the autosampler and the columns (steel, PEEK) or other flushing solutions.
Empty sample contaminated	Use new empty sample.
Cause unknown.	Attempt to solve problem by using different solvents.

If no injection is performed:

Cause of fault	Elimination
Liquid path blocked	<ol style="list-style-type: none"> 1. Disconnect the plastic capillary of the needle from the injection valve. 2. Start system flushing. 3. If solvent escapes at the injection valve connection to the needle, check the needle. 4. If no solvent escapes at the injection valve connection to the needle, disconnect the buffer tube from the injection valve. 5. Start system flushing. 6. If solvent flows out at the open end of the buffer tube, check the rotor seal. 7. If no solvent flows out of the open end of the buffer tube, disconnect the buffer tube from the syringe valve. 8. Start system flushing. 9. If solvent flows out of syringe valve, check the buffer tube. 10. If no solvent flows out of the syringe valve, check whether the connections of the liquid path have been tightened too much.
Leaking valve	<ol style="list-style-type: none"> 1. Unscrew the plastic capillary leading to the needle from the injection valve. 2. Disconnect the plastic capillary leading to the syringe from the injection valve. 3. Connect the HPLC pump to the injection valve. 4. Close the connection to the column at the injection valve. 5. Start the pump at a low flow rate. 6. Check the connections to the syringe and to the needle at the injection valve to ensure that they are tight. 7. If liquid escapes, check the rotor seal. 8. If no liquid escapes, check the HPLC with a manual valve.

System Messages from OpenLAB®

The various system messages of the KNAUER OpenLAB® chromatography software are explained below. The system messages are sorted alphabetically.

System message	Explanation
<i>Autosampler is in run mode.</i>	<ul style="list-style-type: none"> ▪ Quit the control software and restart. ▪ Restart the device.
<i>Autosampler is not responding. Please check communication settings and ensure the device is online.</i>	Restart the device. Check the network settings. Inform the Technical Support of the manufacturer in case the system message repeats itself.
<i>Cannot run autosampler.</i>	Restart the device. Check the network settings. Inform the Technical Support of the manufacturer in case the system message repeats itself.
<i>Cannot set destination vial to %d.</i>	Check parameters in control software and correct entry.
<i>Cannot set first transport vial to %d.</i>	Check parameters in control software and correct entry.
<i>Cannot set last transport vial to %d.</i>	Check parameters in control software and correct entry.
<i>Cannot stop autosampler.</i>	Check the network settings. Inform the Technical Support of the manufacturer in case the system message repeats itself.
<i>Communication port for autosampler was not initialized. Please check the configuration settings.</i>	Restart the device. Inform the Technical Support of the manufacturer in case the system message repeats itself.
<i>Configuration settings do not match with the device. Run cannot start.</i>	Check configuration and settings.
<i>Destination position not reached.</i>	Restart the device. Inform the Technical Support of the manufacturer in case the system message repeats itself.
<i>Deviation of more than +/- 2 mm towards home.</i>	<ul style="list-style-type: none"> ▪ Look for visible obstructions in area of vial plate. ▪ Check the belt tension of the vial plate.

System message	Explanation
<i>Dispenser error.</i>	Restart the device. Inform the Technical Support of the manufacturer in case the system message repeats itself.
<i>Electronics error.</i>	Restart the device. Inform the Technical Support of the manufacturer in case the system message repeats itself.
<i>EEPROM error in adjustments.</i>	Restart the device. Inform the Technical Support of the manufacturer in case the system message repeats itself.
<i>EEPROM error in log counter.</i>	Restart the device. Inform the Technical Support of the manufacturer in case the system message repeats itself.
<i>EEPROM error in settings.</i>	Restart the device. Inform the Technical Support of the manufacturer in case the system message repeats itself.
<i>EEPROM write error.</i>	Restart the device. Inform the Technical Support of the manufacturer in case the system message repeats itself.
Error 369	Not enough transport liquid in store. Refill transport liquid.
Error 370	Not enough reagent in store. Refill reagent.
<i>Error by setting Mix&Dilute vials.</i>	Restart the device. Inform the Technical Support of the manufacturer in case the system message repeats itself.
<i>Error occurred during initialization, the Autosampler AS-1 cannot start.</i>	Restart the device. Inform the Technical Support of the manufacturer in case the system message repeats itself.
<i>Error resetting output.</i>	Restart the device. Inform the Technical Support of the manufacturer in case the system message repeats itself.
<i>Error running user defines program.</i>	Restart the device. Inform the Technical Support of the manufacturer in case the system message repeats itself.

System message	Explanation
<i>Error setting injection mode.</i>	Restart the device. Inform the Technical Support of the manufacturer in case the system message repeats itself.
<i>Error setting needle height.</i>	Restart the device. Inform the Technical Support of the manufacturer in case the system message repeats itself.
<i>Error setting injection mode.</i>	Restart the device. Inform the Technical Support of the manufacturer in case the system message repeats itself.
<i>Error setting syringe speed.</i>	Restart the device. Inform the Technical Support of the manufacturer in case the system message repeats itself.
<i>Error setting the analysis time.</i>	Restart the device. Inform the Technical Support of the manufacturer in case the system message repeats itself.
<i>Error setting the auxiliaries.</i>	Restart the device. Inform the Technical Support of the manufacturer in case the system message repeats itself.
<i>Error setting the flush time.</i>	Restart the device. Inform the Technical Support of the manufacturer in case the system message repeats itself.
<i>Error setting the flush volume.</i>	Restart the device. Inform the Technical Support of the manufacturer in case the system message repeats itself.
<i>Error setting the injection volume.</i>	Restart the device. Inform the Technical Support of the manufacturer in case the system message repeats itself.
<i>Error setting the loop volume.</i>	Restart the device. Inform the Technical Support of the manufacturer in case the system message repeats itself.
<i>Error setting the prep. mode.</i>	Restart the device. Inform the Technical Support of the manufacturer in case the system message repeats itself.

System message	Explanation
<i>Error setting the syringe volume.</i>	Restart the device. Inform the Technical Support of the manufacturer in case the system message repeats itself.
<i>Error setting timed events.</i>	Restart the device. Inform the Technical Support of the manufacturer in case the system message repeats itself.
<i>Error setting the tray configuration.</i>	Restart the device. Inform the Technical Support of the manufacturer in case the system message repeats itself.
<i>Error setting the tray temperature.</i>	Restart the device. Inform the Technical Support of the manufacturer in case the system message repeats itself.
<i>Error setting the vial number.</i>	Restart the device. Inform the Technical Support of the manufacturer in case the system message repeats itself.
<i>Error setting tubing volume.</i>	Restart the device. Inform the Technical Support of the manufacturer in case the system message repeats itself.
<i>Error setting wash volume.</i>	Restart the device. Inform the Technical Support of the manufacturer in case the system message repeats itself.
<i>Flush volume error.</i>	Check parameters in control software and correct entry.
<i>Home sensor activated when not expected.</i>	<ul style="list-style-type: none"> ▪ Check parameters in control software and correct entry. ▪ Restart the device. Inform the Technical Support of the manufacturer in case the system message repeats itself.
<i>Home sensor not de-activated.</i>	<ul style="list-style-type: none"> ▪ Check whether there are visible obstructions impairing the vial plate. ▪ Restart the device. Inform the Technical Support of the manufacturer in case the system message repeats itself.

System message	Explanation
<i>Home sensor not reached.</i>	<ul style="list-style-type: none"> ▪ Check whether there are visible obstructions impairing the vial plate. ▪ Restart the device. Inform the Technical Support of the manufacturer in case the system message repeats itself.
<i>Horizontal: home sensor activated when not expected.</i>	<ul style="list-style-type: none"> ▪ Restart the device. Inform the Technical Support of the manufacturer in case the system message repeats itself.
<i>Horizontal: home sensor not de-activated.</i>	<ul style="list-style-type: none"> ▪ Check whether there are visible obstructions impairing the needle unit. ▪ Restart the device. Inform the Technical Support of the manufacturer in case the system message repeats itself.
<i>Horizontal: home sensor not reached.</i>	<ul style="list-style-type: none"> ▪ Check whether there are visible obstructions impairing the needle unit. ▪ Restart the device. Inform the Technical Support of the manufacturer in case the system message repeats itself.
<i>Horizontal: needle position is unknown.</i>	Initialize the needle unit using the control software.
<i>Illegal sensor readout.</i>	Restart the device. Inform the Technical Support of the manufacturer in case the system message repeats itself.
<i>Incorrect first destination vial.</i>	Check parameters in control software and correct entry.
<i>Injection needle unit error.</i>	<ul style="list-style-type: none"> ▪ Check whether there are visible obstructions impairing the needle unit. ▪ Restart the device. Inform the Technical Support of the manufacturer in case the system message repeats itself.
<i>Injection valve or ISS unit error.</i>	Restart the device. Inform the Technical Support of the manufacturer in case the system message repeats itself.

System message	Explanation
<i>Injection volume %.2f is invalid. For specified injection method, volume should be within the range %.2f µl-%.2f µl, with %.2f. µl increments.</i>	Check parameters in control software and correct entry.
<i>Injection volume error.</i>	Check parameters in control software and correct entry.
<i>Invalid %s vial position %02d. The vial position must be between 01 and %02d.\n.</i>	Check parameters in control software and correct entry.
<i>Invalid combination of the trays. The combination of different trays for the Mix&Dilute mode is not allowed.</i>	<ul style="list-style-type: none"> ▪ Insert the correct vial plate. ▪ Check parameters in control software and correct entry.
<i>Invalid combination of the trays. The combination of plates 384 low and 96 high is not allowed.</i>	<ul style="list-style-type: none"> ▪ Insert the correct vial plate. ▪ Check parameters in control software and correct entry.
<i>Invalid configuration. ISS option not installed on autosampler. Please switch off this option in configuration dialog.</i>	Check parameters in control software and correct entry.
<i>Invalid configuration. SSV option not installed on autosampler. Please switch off this option in configuration dialog.</i>	Check parameters in control software and correct entry.
<i>Invalid flush volume %d µl. The flush volume should be between 0 and %d µl.</i>	Check parameters in control software and correct entry.
<i>Invalid flush volume %.2f µl. The flush volume should be between 0 and %.2f µl.</i>	Check parameters in control software and correct entry.
<i>Invalid input. Only values with increments of %.2f allowed.</i>	Check parameters in control software and correct entry.
<i>Invalid integer number.</i>	Check parameters in control software and correct entry.

System message	Explanation
<i>Invalid instrument is detected.</i>	Check parameters in control software and correct entry.
<i>Invalid loop volume %d µl. The loop volume should be between 0 and %d µl.</i>	Check parameters in control software and correct entry.
<i>Invalid loop volume %2f µl. The loop volume should be between 0 and %2f µl.</i>	Check parameters in control software and correct entry.
<i>Invalid mix program: no Destination vial is specified in the configuration dialog.</i>	Check parameters in control software and correct entry.
<i>Invalid mix program: no Reagent A vial is specified in the configuration dialog.</i>	Check parameters in control software and correct entry.
<i>Invalid mix program: no Reagent B vial is specified in the configuration dialog.</i>	Check parameters in control software and correct entry.
<i>Invalid mix times. The time should be between 1 and 9.</i>	Check parameters in control software and correct entry.
<i>Invalid needle height %d mm. The needle height should be between %d and %d mm.</i>	Check parameters in control software and correct entry.
<i>Invalid time-based method. Several AUX events have the same time.</i>	Check parameters in control software and correct entry.
<i>Invalid time-based method. Several SSV events have the same time.</i>	Check parameters in control software and correct entry.
<i>Invalid tray temperature %d °C. The temperature should be between 4 and 22 °C.</i>	Check parameters in control software and correct entry.
<i>Invalid tray configuration: two or more vial positions are the same.</i>	Check parameters in control software and correct entry.

System message	Explanation
<i>Invalid tubing volume %d µl. The tubing volume should be between 0 and %d µl.</i>	Check parameters in control software and correct entry.
<i>Invalid tubing volume %2f µl. The tubing volume should be between %2f and %2f µl.</i>	Check parameters in control software and correct entry.
<i>Invalid wait time. The time should be between 0 and 9 h 50 min 59 sec. Invalid wash volume %d µl. The wash volume should be between %d and %d µl</i>	Check parameters in control software and correct entry.
<i>Invalid volume %d µl. The volume should be between the 0 and the syringe volume (%d µl).</i>	Check parameters in control software and correct entry.
<i>ISS valve error.</i>	Restart the device. Inform the Technical Support of the manufacturer in case the system message repeats itself.
<i>Missing destination vial.</i>	<ul style="list-style-type: none"> ▪ Check position of sample vial. ▪ Check parameters in control software and correct entry.
<i>Missing reagent vial.</i>	<ul style="list-style-type: none"> ▪ Check position of sample vial. ▪ Check parameters in control software and correct entry.
<i>Missing transport vial.</i>	<ul style="list-style-type: none"> ▪ Check position of sample vial. ▪ Check parameters in control software and correct entry.
<i>Needle movement error.</i>	<ul style="list-style-type: none"> ▪ Check position of needle unit. ▪ Restart the device.
<i>Missing vial.</i>	<ul style="list-style-type: none"> ▪ Check position of needle unit. ▪ Restart the device.
<i>Missing wash vial error.</i>	<ul style="list-style-type: none"> ▪ Check position of needle unit. ▪ Restart the device.
<i>No destination vial is specified in the configuration.</i>	Check parameters in control software and correct entry.
<i>No reagent A vial is specified in the configuration.</i>	Check parameters in control software and correct entry.

System message	Explanation
<i>No reagent B vial is specified in the configuration.</i>	Check parameters in control software and correct entry.
<i>No transport vials are defined in the tray configuration. It is not possible to use the μl pick-up injection mode.</i>	Check parameters in control software and correct entry.
<i>No user defined or mix program is running.</i>	Check parameters in control software and correct entry.
<i>Not enough reagent liquid.</i>	Check volume of liquid and change as required.
<i>Not enough transport liquid available due to missing transport vials.</i>	Check volume of liquid and change as required.
<i>Please specify inject marker or AUX event to be able to trigger the run.</i>	Check parameters in control software and correct entry.
<i>Selecting transport position failed.</i>	Restart the device. Inform the Technical Support of the manufacturer in case the system message repeats itself.
<i>Serial number is not valid. Please check the configuration.</i>	Check parameters in control software and correct entry.
<i>Setting mix program error.</i>	Restart the device. Inform the Technical Support of the manufacturer in case the system message repeats itself.
<i>Setting service mode failed.</i>	Restart the device. Inform the Technical Support of the manufacturer in case the system message repeats itself.
<i>Syringe dispenser unit error.</i>	Restart the device. Inform the Technical Support of the manufacturer in case the system message repeats itself.
<i>Syringe home sensor not de-activated.</i>	<ul style="list-style-type: none"> ▪ Needle flushing with control software. ▪ Restart the device. Inform the Technical Support of the manufacturer in case the system message repeats itself.

System message	Explanation
<i>Syringe home sensor not reached.</i>	<ul style="list-style-type: none"> ▪ Needle flushing with control software. ▪ Restart the device. Inform the Technical Support of the manufacturer in case the system message repeats itself.
<i>Syringe position is unknown.</i>	Initialize the syringe unit using the control software.
<i>Syringe rotation error.</i>	<ul style="list-style-type: none"> ▪ Needle flushing with control software. ▪ Restart the device. Inform the Technical Support of the manufacturer in case the system message repeats itself.
<i>Syringe valve did not find destination position.</i>	<ul style="list-style-type: none"> ▪ Needle flushing with control software. ▪ Restart the device. Inform the Technical Support of the manufacturer in case the system message repeats itself.
<i>Temperature above 48 °C at cooling ON.</i>	<ul style="list-style-type: none"> ▪ Switch off the cooling and check whether ambient temperature sensor is properly functioning. ▪ Inform the Technical Support of the manufacturer in case the system message repeats itself.
<i>The 10 ml syringe cannot be used for standard injections.</i>	Exchange the syringe.
<i>The autosampler has detected another tray than that which is currently configured. Please select the correct tray in configuration dialog.</i>	Check control software configuration and correct entry.
<i>ISS-A option not installed on autosampler. Please switch off ISS-A option in configuration dialog.</i>	Check control software configuration and correct entry.
<i>ISS-B option not installed on autosampler. Please switch off ISS-B option in configuration dialog.</i>	Check control software configuration and correct entry.

System message	Explanation
<i>Oven option not installed on autosampler. Please switch off oven option in configuration dialog.</i>	Check control software configuration and correct entry.
<i>The autosampler is not ready. Please try later.</i>	Restart the device. Inform the Technical Support of the manufacturer in case the system message repeats itself.
<i>The injection volume of %2f µl is invalid. For the specified injection method, volume should equal %2f µl.</i>	Check parameters in control software and correct entry.
<i>The sample needle is not in the home position while the tray is rotating.</i>	Check parameters in control software and correct entry.
<i>Trace from tray cooling cannot be acquired. Tray cooling is off.</i>	Check parameters in control software and correct entry.
<i>Tray advance is not available at this time.</i>	Check parameters in control software and correct entry.
<i>Tray error.</i>	Check parameters in control software and correct entry.
<i>Valve error.</i>	Check parameters in control software and correct entry.
<i>Vertical: home sensor not de-activated.</i>	<ul style="list-style-type: none"> ▪ Check whether there are visible obstructions impairing the needle unit. ▪ Restart the device. Inform the Technical Support of the manufacturer in case the system message repeats itself.
<i>Vertical: home sensor not reached.</i>	<ul style="list-style-type: none"> ▪ Check whether there are visible obstructions impairing the needle unit. ▪ Restart the device. Inform the Technical Support of the manufacturer in case the system message repeats itself.
<i>Vertical: needle position is unknown.</i>	Initialize the instrument in the control software.

System message	Explanation
<i>Vertical: stripper did not detect plate (or wash/waste). Missing vial.</i>	<ul style="list-style-type: none"> ▪ Check sample vial and plate. ▪ Restart the device. Inform the Technical Support of the manufacturer in case the system message repeats itself.
<i>Vertical: stripper stuck.</i>	Restart the device. Inform the Technical Support of the manufacturer in case the system message repeats itself.
<i>Vertical: The sample needle arm is at an invalid position.</i>	Restart the device. Inform the Technical Support of the manufacturer in case the system message repeats itself.
<i>Vial number error.</i>	Check parameters in control software and correct entry.
<i>Wear-out limit reached.</i>	Restart the device. If the system message appears again, notify KNAUER Technical Support. The valve must be replaced.
<i>Wrong loop volume. The largest loop volume for standard injections is 1000 µl.</i>	Check parameters in control software and correct entry.
<i>Wrong tubing volume. The largest tubing volume for standard injections is 200 µl.</i>	Check parameters in control software and correct entry.

Technical Data

Ambient conditions	Temperature range	10 – 40 °C; 50 – 104 °F
	Humidity	20–80 % RH

Standard Version

Autosampler 3950	Injection pressure	Up to 1000 bar
	Sample capacity	Microtiter plates for max. 768 wells or sample plates for max. 96 vials
	Injection volume	0.1-5000 µl
	Sample loop	10 µl
	Syringe	250 µl, 500 µl
	Injection time	<60 s, incl. cleaning
	Injection modes	PASA™ loop injection principle: <ul style="list-style-type: none"> ▪ Full loop injection ▪ Partial loop injection ▪ <i>Microliter pick-up</i> injection
	Reproducibility	RSD (Relative Standard Deviation): <ul style="list-style-type: none"> ▪ Full loop injection < 0.3 % ▪ Partial loop fill with injection volumes > 10 µl: < 0.5 % ▪ Microliter pickup injections with injection volumes > 10 µl: < 1.0 %
	Carryover	< 0.05 % with needle cleaning
	Temperature control	4 – 40 °C ¹
	Weight	<ul style="list-style-type: none"> ▪ 19 kg ▪ With temperature control: 21 kg
	Dimensions (length x width x height)	<ul style="list-style-type: none"> ▪ 510 × 300 × 360 mm ▪ With cooling: 575 × 300 × 360 mm
	Supply voltage range	95-240 V
	Supply frequency	50-60 Hz
Pollution degree	2	
Installation category	II	

1. Relevant to Autosampler 3950 with temperature control

Product Range

Device and Accessories

Name	Order number
Autosampler 3950, LAN, Standard	A50070
Autosampler 3950, LAN, cool/heat	A500701
User manual	V1511

Repeat Orders

Name	Order number
Accessory kit for autosampler with vials, caps, septum and pliers for opening and closing	A0664
Vial plates for 1.5 ml vials 2 pieces	A50050
96 well plate, U sanitized, 0.35 ml	A1823
96 well plate, U sanitized, 1.2 ml	A1823V1
PTFE tubing, 3.2 mm (1/8") OD, 1.5 mm ID, 300 cm length	A0732
Wash bottle 250 ml	M2054
Syringe 500 µl	M2070
Air needle, white, 62 mm	A50058
Set of air needles	A50059
Air needle, yellow, 50 mm	M20401
Air needle, red, 56 mm	M20402
Air needle, blue, 68 mm	M20403
Air needle, green, 74 mm	M20404
Air needle, black, 80 mm	M20405
2 x 2.5 A fuse	M2067
Network cable	A5255
RS-232 cable for connection to a PC	A0895
User manual, German	V1511
User manual, English	V1511A

Accessories 84 +3

Name	Order number
Vial plate for 84x1.5ml and 3x10ml vials	A500501
125 Vials 10 ml, 500 crimp caps and 500 septa, ø 22 mm	A1662

Legal Information

Warranty Conditions

The factory warranty for the device is valid for 12 months after the date of dispatch. All warranty claims shall expire in the event that any unauthorized changes are made to the device.

During the warranty period, any components with material or design-related defects will be replaced or repaired by the manufacturer free of charge.

This warranty excludes the following:

- Accidental or willful damage
- Damage or errors caused by third parties that are not contractually related to the manufacturer at the time the damage occurs
- Wear parts, fuses, glass parts, columns, light sources, cuvettes and other optical components
- Damage caused by negligence or improper operation of the device and damage caused by clogged capillary
- Packaging and transport damage

In the event of device malfunctions, directly contact the manufacturer.

Wissenschaftliche Gerätebau
Dr. Ing. Herbert Knauer GmbH
Hegauer Weg 38
14163 Berlin, Germany
Phone: +49 30 809727-111
Telefax: +49 30 8015010
e-mail: info@knauer.net
Internet: www.knauer.net

Transport Damage

The packaging of our devices provides the best possible protection against transport damage. Check the devices for signs of transport damage. In case you notice any damage, contact the Technical Support and the forwarder company within three workdays.

HPLC Glossary

In the following chapter you find abbreviations and terminology that is used in HPLC.

Term	Definition
analytical	The analysis and quantitative measuring of samples in HPLC (see preparative)
capillary	Thin metal or PEEK pipe that connects the components and devices inside of an HPLC system.
chromatogram	The record of a detector signal, depending on the flow volume of the mobile phase or time.
column	Pipe with caps which are permeated by the solvent. The pipe contains the packing material.
detector	Device measuring the composition or the quantity of a substance.
GLP	Quality assurance system for laboratories (Good Laboratory Practice)
HPLC	High-pressure liquid chromatography (HPLC). High-pressure liquid chromatography
pump	Device which conveys solvents in controlled quantities into a system.
sample	A mixture of different components which are to be separated by chromatography. They are transported by the mobile phase and removed from the column.
sample loop	A loop that is separated from the chromatographic system through a valve and which receives the sample first. By switching the valve, the stream of solvent is lead through the loop and the sample is flushed onto the separating column.
solvent	Flowing agent that transports the substances to be separated or isolated through the column (solvent, mobile phase)
valve	The mechanism which the sample is injected into the solvent stream.

List of Figures

Fig. 1:	Autosampler 3950, with optional temperature control and coolunit	7
Fig. 2:	System adapter for autosampler	18
Fig. 3:	System adapter installed	18
Fig. 4:	Front view of the device	19
Fig. 5:	Push the door into the interior	20
Fig. 6:	Removing the front panel	20
Fig. 7:	Removing the coolunit cover	20
Fig. 8:	Front view with sample compartment	21
Fig. 9:	Rear panel of the device	22
Fig. 10:	ILD™ injection principle	25
Fig. 11:	Full loop injection: Initial conditions	27
Fig. 12:	Full loop injection: The needle and sample lines are flushed	27
Fig. 13:	Full loop injection: Valve switches to <i>LOAD</i> position	28
Fig. 14:	Full loop injection: The sample loop is fully filled.	28
Fig. 15:	Full loop injection: Valve switches to <i>INJECT</i> position	29
Fig. 16:	Air Segment for full loop injection	29
Fig. 17:	Air segment for partial loop full	30
Fig. 18:	Scheme of initial position for partial loop fill	30
Fig. 19:	Partial loop fill: Valve position "Inject"	31
Fig. 20:	Partial loop fill: Aspirating the sample	31
Fig. 21:	Partial loop fill: Valve position "Load"	31
Fig. 22:	Partial loop fill: Injecting into the column	32
Fig. 23:	Partial loop fill: Aspirating flushing solution	32
Fig. 24:	Partial loop fill: Flushing sample and solvent	32
Fig. 25:	Microliter pick-up injection: Initial conditions	33
Fig. 26:	Microliter pick-up injection: Sample line is filled with solvent	33
Fig. 27:	Microliter pick-up: Sample is aspirated	34
Fig. 28:	Microliter pick-up: The injection volume is aspirated.	34
Fig. 29:	Microliter pick-up: The sample is transported through the sample loop.	35
Fig. 30:	Microliter pick-up: The sample is transported to the column.	35
Fig. 31:	Microliter pick-up injection with air segment (A), without air segment (B)	36
Fig. 32:	Microliter pick-up 84+3: Initial conditions	36
Fig. 33:	Microliter pick-up 84+3: Sample line is filled with flushing solution.	36
Fig. 34:	Microliter pick-up 84+3: Injection valve switches to <i>LOAD</i> position	37
Fig. 35:	Microliter pick-up 84+3: The sample is transported through the sample loop.	37
Fig. 36:	Microliter pick-up 84+3: Injection valve switches to position <i>INJECT</i> .	37
Fig. 37:	Microliter pick-up injection 84+3 with air segment (A),	

	without air segment (B)	38
Fig. 38:	Height of 84+3 vial plate	39
Fig. 39:	Needle length and liquid level	39
Fig. 40:	Sequence	39
Fig. 41:	Reagent/Transport position on 84+3 vial plate	40
Fig. 42:	Standard air needle with 10 and 2 ml sample vial	42
Fig. 43:	Standard air needle with microtiter plates	42
Fig. 44:	Calculating the correct air needle	43
Fig. 45:	Selected air needles with sample vials	44
Fig. 46:	Selected air needles with microtiter plates	44
Fig. 47:	Sample positions in columns	47
Fig. 48:	Sample positions in rows	48
Fig. 49:	Reagent/Transport position on 84+3 vial plate	49
Fig. 50:	Connection scheme for capillaries and tubing	50
Fig. 51:	Connections on the injection valve	51
Fig. 52:	Syringe connections	51
Fig. 53:	Tube guide for flushing solution tube	52
Fig. 54:	Connecting the drainage tubing	53
Fig. 55:	Autosampler configuration with ClarityChrom®	54
Fig. 56:	System flushing with the Service Manager software	56
Fig. 57:	Operation Qualification Report	65
Fig. 58:	Disassembling the valve	68
Fig. 59:	Removing the rotor seal	68
Fig. 60:	Rotor and rotor seal	68
Fig. 61:	Installing the valve	69
Fig. 62:	Exchanging the sample needle	71
Fig. 63:	Exchanging the air needle	72
Fig. 64:	Exchanging the syringe	73
Fig. 65:	Exchanging the syringe valve	75

Index

A

- Accessories 15
- Accessories kit 15
- Additives 9
- Air needle 26, 42
 - calculation example 45
 - choosing the correct air needle 43
 - standard air needle 42
 - types 44
- Air segment
 - full loop injection 29
- analytical 98
- Auto-injection system 19
- Automatic configuration 23
- AVV marking 77

C

- Capillary 98
- Capillary connector 50
- Care 66, 76
- Chromatogram 98
- Cleaning 76
- Column 98
- Connecting the autosampler to other devices 54
- Connecting the valve 51
- Connection
 - closed-contact output 57
 - drainage tubing 53
 - I/O connection 56, 57
 - LAN 54
 - syringe 51
 - TTL inputs 56
 - valve 51
- Connection to computer 23
- Coolunit 7

D

- Decontamination 13
- Degasser 10
- Detector 98
- Device Overview 7

- Device test 58
 - archiving 63
 - linearity 61
 - mixture test 62
 - reproducibility 59
 - sample carryover 60
 - test interval 58

E

- Environmental protection 77

F

- Features 7
- Front view of the device 19, 21
- Full loop injection 25, 26, 27
 - air segment 29

G

- GLP 98

H

- Handling 45
 - sample vial 45
- Hose connectors 50
- HPLC 98
- HPLC Glossary 98

I

- ILD™ 25
- Initial Startup 19
- Injection methods 25
- Injection principle
 - ILD™ 25
 - loop injection 25
- Injection principles 25
- Installation 15
- Installation location, see location 16
- Installing a local area network 23
- Intended Use 7

L

- LAN 23
 - connection problems 78
 - connectivity problems 23
- Leak 12
- Legal Information 97
- List of figures 99
- Location 16
- Loop fill
 - partial 30

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Function Test

No.	Test	Settings	Specification	Result
1	Reproducibility	For the sample volume reproducibility the prepared uracil solution a (50 ppm) is injected again with an injection volume of 2 µl.	The variation coefficient may not exceed 2 %.	Test successful 1,229 %
2	Carryover	Inject Uracil and Eluent by turns. Calculate sample carry-over percent.	The sample carryover percent may not exceed 1 %.	Test successful 0,944 %
3	Linearity	Injection of 5 test solutions.	Correlation coefficient R ≥ 0,97	Test successful R = 1,000

Date of inspection: 6 January 2010

Next inspection on: _____

Inspector: _____

Signature: _____

Abb. 58 Operation Qualification Report

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All rights reserved.
® AZURA is a registered trademark of