

NEWS



Individually dosed –
Therapeutic drug monitoring with LC/MS/MS

Bread packaging –
What's hidden in plastic bread packaging?

Hair in your soup?
Microchip electrophoresis can help

Acid Test – TOC determination in concentrated hydrochloric acid

Headspace-Cold Trap Sampling




Figure 1: Left: Peak of $m/z = 62$ (vinyl chloride) for different cryofocus temperatures (without cryofocus, $-20\text{ }^{\circ}\text{C}$, $-70\text{ }^{\circ}\text{C}$, $-130\text{ }^{\circ}\text{C}$ and $-140\text{ }^{\circ}\text{C}$). Right: Peak of $m/z = 62$ for different heating rates of the cryofocus after refocusing.

Hans-Ulrich Baier, Panos Meletis und Stephan Schröder, Shimadzu Deutschland, Duisburg, Germany

Analysis of EPA624 regulated volatile organic compounds in drinking and wastewater is usually performed with headspace or purge and trap technique using a so-called 624 phase with 30 m, 0.25 mm and 1.4 μm . Reducing analysis time (fast GC) but maintaining chromatographic resolution has been applied successfully using narrow bore columns in various fields. However, the reported results were based mainly on liquid injection techniques.

In headspace analysis the transfer of sample from the insert to the column is quite slow as small split ratios are normally used in favor of sensitivity. The spatial distribution of analyte molecules in the glass insert cannot therefore be refocused easily, and fast GC approaches are difficult. A cold trap (cryofocus, ATASGL,

The Netherlands) was therefore mounted at the top of the column directly under the injector, cooling the first part of the column in order to refocus volatile compounds showing a broad band during passage through the injector liner. The cooling was established by direct transfer of liquid nitrogen to the trap.

Instrumentation used was a Shimadzu GCMS-QP2010 Ultra with an AOC-5000 Plus headspace sampler. As the column is surrounded by the directly cooled cryofocus, refocusing takes place

inside the column. In this study, the inner diameter of the chosen column was 0.18 mm. Length and film thickness were 20 m and 1 μm respectively. Split ratio was 5:1 and the linear velocity was




Figure 2: Full scan chromatogram (TIC) of 60 volatile compounds

Sampling

Fast GCMS analysis of VOCs in water




Figure 3: Calibration curves for benzene and vinyl chloride

set to 45 cm/s. The GC oven temperature began at 40 °C for 5 min and then ramped with 50 °C/min to 120 °C, 30 °C/min to 170 °C, 60 °C/min to 220 °C. Injection volume was 1 mL headspace from a 20 mL vial filled with 5 mL water matrix. Different cold trap temperatures were set. The mass spectrometer was operated in scan and selected ion monitoring (SIM) mode for highly sensitive analysis.

Good peak shapes at 50 °C/s

Figure 1 (left) shows the m/z 62 relative to vinyl chloride for different cryofocus temperatures. The largest effect of refocusing monitored by measuring the peak profile at the end of the column in the mass spectrometric detector was observed at -140 °C cold trap temperature and subsequent heating to 250 °C at a rate of 50 °C/s. Figure 1 (right) shows the influence of different heating rates for vinyl chloride indicating that 50 °C/s ensures that the re-

leasing process is fast enough to obtain good peak shapes. The peakwidth at half maximum is 8 sec and 0.5 sec for cryofocus temperatures of 0 °C and -140 °C, respectively. The peak height is increased drastically, with considerable improvement to the limit of detection (LOD). The complete chromatogram is shown in figure 2 and the compounds are listed in table 1.

Analysis within ten minutes

Analysis time for 60 volatile compounds was less than ten minutes. Calibration was performed between 0.001 and 1 µg/L. The regression coefficient R showed values of better than 0.998 for all compounds, indicating the high method precision. Two curves are shown in figure 3.

The LOD for benzene and vinyl chloride was determined as 0.005 µg/L and 0.001 µg/L respectively. In figure 4, the selected ion mass traces for tetrachloroethene and




Figure 4: Peaks of tetrachloroethene and 1,1,2-trichloroethene measured for a water sample taken from the Rhine River

ID#	Name	Ret. Time			
1	Dichlorodifluoromethane	0.945	30	Dibromochloromethane	7.260
2	Chloromethane	1.077	31	1,2-Dibromoethane	7.319
3	Vinyl chloride	1.156	32	Chlorobenzene	7.547
4	Bromomethane	1.375	33	Ethylbenzene	7.582
5	Chloroethane	1.443	34	1,1,1,2-Tetrachloroethane	7.586
6	Trichlorofluoromethane	1.587	35	p-Xylene	7.639
7	1,1-Dichloroethene	1.970	36	m-Xylene	7.639
8	Methylene chloride	2.432	37	o-Xylene	7.835
9	trans-1,2-Dichloroethene	2.645	38	Styrene	7.849
10	1,1-Dichloroethane	3.163	39	Tribromomethane	7.957
11	2,2-Dichloropropane	3.890	40	Isopropylbenzene	8.002
12	cis-1,2-Dichloroethene	3.967	41	Bromobenzene	8.177
13	Bromochloromethane	4.346	42	1,1,2,2-Tetrachloroethane	8.186
14	Trichloromethane	4.521	43	1,2,3-Trichloropropane	8.216
15	1,1,1-Trichloroethane	4.683	44	n-Propylbenzene	8.200
16	Tetrachloromethane	4.900	45	2-Chlorotoluene	8.262
17	1,1-Dichloropropene	5.005	46	1,3,5-Trimethylbenzene	8.286
18	Benzene	5.313	47	4-Chlorotoluene	8.319
19	1,2-Dichloroethane	5.475	48	tert-Butylbenzene	8.438
20	Trichloroethene	6.000	49	1,2,4-Trimethylbenzene	8.466
21	1,2-Dichloropropane	6.197	50	sec-Butylbenzene	8.538
22	Dibromomethane	6.284	51	4-Isopropyltoluene	8.599
23	Bromodichloromethane	6.414	52	1,3-Dichlorobenzene	8.616
24	cis-1,3-Dichloropropene	6.675	53	1,4-Dichlorobenzene	8.658
25	Toluene	6.813	54	n-Butylbenzene	8.780
26	trans-1,3-Dichloropropene	6.992	55	1,2-Dichlorobenzene	8.825
27	Tetrachloroethene	7.084	56	1,2-Dibromo-3-chloropropane	9.159
28	1,1,2-Trichloroethane	7.084	57	1,2,4-Trichlorobenzene	9.489
29	1,3-Dichloropropane	7.165	58	1,1,2,3,4,4-Hexachloro-1,3-b	9.532
30	Dibromochloromethane	7.260	59	Naphthalene	9.612
			60	1,2,3-trichlorobenzene	9.724

Table 1: List of volatile compounds and retention times in minutes

1,1,2-trichloroethane of a real sample (the Rhine River water) is shown. Both concentrations were determined as 0.02 µg/L.