

Comprehensive Listings of Polypropylene
96-Well μLplates® & Piercible Mats

MicroLiter has a product line of polypropylene microplates and piercible covers that offers secure sample containment throughout the analytical sample run. Whether you want to simply stabilize the sample until it's injected or if you want some measure of stability, even after the sample run has finished, we have the solution. To do so we offer two mats, an EVA mat that stabilizes the sample yet is piercible by most instruments and a silicone mat that has the resealing qualities of standard autosampler vial septa.

MicroLiter offers microplates in two configurations, round well and square well. Any round well mat will fit any round well microplate. Likewise, any square well mat will fit any square well microplate.

μLplate® Covers for 96-Well Microplates


Cat. No.	Description	Qty
07-0010N	Cover for Round Well μLplates, EVA	1
07-0011N	Piercible Cover for All Round Well μLplates, EVA	1
07-0015N	Piercible Cover for All Square Well μLplates, EVA	1
07-0017N	Piercible Cover for All Square Well μLplates, EVA	1
07-0060N	Silicone Mat w/ Sprayed-On PTFE Barrier Square Well	1
07-0061N	Silicone Mat w/ Sprayed-On PTFE Barrier, Square Well, Pre-Cut	1
07-0061Y	Yellow Silicone Mat w/ Sprayed-On PTFE Barrier, Square Well, Pre-Cut	1
07-0063N	Silicone Mat w/ Sprayed-On PTFE Barrier, Round Well, Pre-Cut	1
07-0065N	Silicone Mat w/ Sprayed-On PTFE Barrier, Round Well	1
07-0066N	Silicone Mat w/ Sprayed-On PTFE Barrier, Round Well, Pre-Cut, Thin	1
07-0066N-B	Silicone Mat w/ Sprayed-On PTFE Barrier, Round Well, Pre-Cut, Bulk	10
07-0067N	Silicone Mat w/ Sprayed-On PTFE Barrier, Square Well, Thin	1

Shallow & Medium μLplates®, Polypropylene


Cat. No.	Description	Qty
07-1000	350μL, Polypropylene, Shallow, Round Well / Flat Bottom	1
07-1100	350μL, Polypropylene, Shallow, Round Well / Round Bottom	1
07-1200	350μL, Polypropylene, Shallow, Round Well / Conical Bottom	1
07-2000	500μL, Polypropylene, Medium, Round Well / Round Bottom	1
07-2100	500μL, Polypropylene, Medium, Round Well / Round Bottom	1
07-2100-B	500μL, PP, Medium, Round Well / Round Bottom, Bulk	10

Deep μLplates®, Polypropylene


Cat. No.	Description	Qty
07-3000	1.3mL, Polypropylene, Deep, Round Well / Round Bottom	1
07-3400	1.3mL, Polypropylene, Deep, Round Well / Round Bottom	1
07-5000	350μL, Polypropylene, Deep, Square Well / Conical Bottom	1
07-6000	1.0mL, Polypropylene, Deep, Square Well / Conical Bottom	1
07-7000	2.0mL, Polypropylene, Deep, Square Well / Conical Bottom	1
07-7400	2.0mL, Polypropylene, Deep, Square Well / Conical Bottom	1
07-8000	2.0mL, Polypropylene, Deep, Square Well / Round Bottom	1
07-8400	2.0mL, Polypropylene, Deep, Square Well / Round Bottom	1




07-0010N




07-0011N




07-0015N & 07-0017N




07-0060N




07-0063N




07-0065N




07-0066N



07-0067N




07-1000




07-1100




07-1200




07-2000 & 07-2100




07-3000




07-3400




07-5000



07-6000



07-7000 & 07-7400



07-8000 & 07-8400