







Chart of Pipette Tips






Microtips, 10µL

Reference	Cat No.	Length (mm)	Description	Qty / Case	DNase, RNase Free	Filter	Pyrogen Free	Sterile
 A	851172	45	Short Tip, 10 racks of 100, Natural	1000	✓			
 B	851174	76	Long Tip, 4 racks of 200, Natural	800	✓			
 E	W851239-01	31	10 racks of 96 tips, Natural	960				
 F	W851239-02	31	10 racks of 96 tips, Natural	960	✓		✓	✓
 G	W851239-03	31	Bulk pack of 1000 tips, Natural	1000				





Microtips, 100µL

Reference	Cat No.	Length (mm)	Description	Qty / Case	DNase, RNase Free	Filter	Pyrogen Free	Sterile
 H	851179-02	54	10 racks of 96 tips, Natural	960	✓	✓	✓	✓





Microtips, 200µL

Reference	Cat No.	Length (mm)	Description	Qty / Case	DNase, RNase Free	Filter	Pyrogen Free	Sterile
 I	851246	50	10 racks of 96 tips, Natural	960	✓			
 J	851247	50	10 racks of 96 tips, Natural	960	✓		✓	✓
 K	851248	50	Bulk pack of 1000 tips, Natural	1000				
 L	851271	53	Bulk pack of 1000 tips, Yellow	1000	✓		✓	
 M	851272	53	10 racks of 96 tips, Yellow	960	✓		✓	





Microtips, 300 / 350µL

Reference	Cat No.	Length (mm)	Description	Qty / Case	DNase, RNase Free	Filter	Pyrogen Free	Sterile
 N	851181-01	59	300µL - Bulk pack of 1000 tips, Natural	1000	✓	✓		
 O	851181-04	59	300µL - 10 racks of 96 tips, Natural	960	✓	✓	✓	✓
 P	W851189	59	350µL - Bulk pack of 1000 tips, Natural	1000				
 Q	W851190	59	350µL - 10 racks of 96 tips, Natural	960				

Microtips, 1000µL

Reference	Cat No.	Length (mm)	Description	Qty / Case	DNase, RNase Free	Filter	Pyrogen Free	Sterile
 R	851276	78	Bulk pack of 1000 tips, Blue	1000	✓		✓	
 S	851277	78	10 racks of 100 tips, Blue	1000	✓		✓	
 T	851180-01	88	Individually wrapped, Natural	1000	✓	✓	✓	✓
 U	851180-02	88	10 racks of 100 tips, Natural	1000	✓	✓	✓	✓

Macrotips, 2, 5, 10mL

Reference	Cat No.	Length (mm)	Description	Qty / Case
 V	851355	115	2mL - Bulk pack of 250 tips, Natural	250
 W	851358	150	10mL - Bulk pack of 100 tips, Natural	100
	W835150	—	Autoclaving Box for 24 - 10mL tips	4
 X ^{1,2}	W851362	123	5mL - 2 packs of 250 tips, Natural*	500
 Y ^{1,2}	W851363	123	5mL - 2 racks of 50 tips, Natural*	100

Polypropylene Straw Tips

Reference	Cat No.	Length (mm)	Description	Qty / Case
AA	851344	190	1.1mL - 4 packs of 25 tips, Natural	250

Cat. No.	Model	Volume	Tip Ref.
851160	822	0.2-2µL	A,B,E,F,G
851162	822	1-10µL	A,B,E,F,G
851163	822	2-20µL	I,J,K,L,M,N,O
851164	822	10-100µL	H,I,J,K,L,M,N,O
851166	822	20-200µL	I,J,K,L,M,N,O
851168	822	100-1000µL	R,S,T,U
851183-2	852	10-100µL	H,I,J,K,L,M,N,O
851184-2	852	10-100µL	H,I,J,K,L,M,N,O
851187-2	852	20-200µL	I,J,K,L,M,N,O
851188-2	852	20-200µL	I,J,K,L,M,N,O
851206	815	10YµL	I,J,K,L,M,N,O
851340	832	0.2-2mL	V
851343	810	0.1 & 1mL	AA
851345	832	1-10mL	W
W810000	815	1µL	A,B,E,F,G
W810002	815	5µL	A,B,E,F,G
W810004	815	10µL	A,B,E,F,G
W810006	815	15µL	I,J,K,L,M,N,O
W810008	815	20µL	I,J,K,L,M,N,O
W810010	815	25µL	H,I,J,K,L,M,N,O
W810012	815	30µL	H,I,J,K,L,M,N,O
W810013	815	32µL	H,I,J,K,L,M,N,O
W810014	815	40µL	H,I,J,K,L,M,N,O
W810016	815	50µL	H,I,J,K,L,M,N,O
W810018	815	60µL	H,I,J,K,L,M,N,O
W810020	815	70µL	H,I,J,K,L,M,N,O
W810022	815	75µL	H,I,J,K,L,M,N,O
W810024	815	80µL	H,I,J,K,L,M,N,O
W810026	815	90µL	H,I,J,K,L,M,N,O
W810028	815	100µL	H,I,J,K,L,M,N,O
W810030	815	120µL	I,J,K,L,M,N,O
W810032	815	150µL	I,J,K,L,M,N,O
W810034	815	200µL	I,J,K,L,M,N,O
W810036	815	250µL	R,S,T,U
W810038	815	300µL	R,S,T,U
W810040	815	400µL	R,S,T,U
W810042	815	500µL	R,S,T,U
W810044	815	600µL	R,S,T,U
W810046	815	700µL	R,S,T,U
W810048	815	750µL	R,S,T,U
W810050	815	800µL	R,S,T,U
W810052	815	900µL	R,S,T,U
W810054	815	1000µL	R,S,T,U
W810056	835F	2mL	V,X,Y
W810058	835F	2.5mL	X,Y
W810060	835F	5mL	X,Y
W810062	835F	10mL	W

Cat. No.	Model	Volume	Tip Ref.
W810300	825	0.2-2µL	A,B,E,F,G
W810300-XS	826	0.2-2µL	A,B,E,F,G
W810302	825	0.5-10µL	A,B,E,F,G
W810302-XS	826	0.5-10µL	A,B,E,F,G
W810304	825	1-10YµL	I,J,K,L,M,N,O
W810304-XS	826	1-10YµL	I,J,K,L,M,N,O
W810306	825	2-20µL	I,J,K,L,M,N,O
W810306-XS	826	2-20µL	I,J,K,L,M,N,O
W810308	825	5-50µL	H,I,J,K,L,M,N,O
W810308-XS	826	5-50µL	H,I,J,K,L,M,N,O
W810310	825	10-100µL	H,I,J,K,L,M,N,O
W810310-XS	826	10-100µL	H,I,J,K,L,M,N,O
W810312	825	20-200µL	I,J,K,L,M,N,O
W810312-XS	826	20-200µL	I,J,K,L,M,N,O
W810314	825	100-1000µL	R,S,T,U
W810314-XS	826	100-1000µL	R,S,T,U
W810315	835	0.2-2mL	V,X,Y
W810316	835	0.5-5mL	X,Y
W810318	835	1-10mL	W
W810600	855	0.5-10µL	A,B,E,F,G
W810602	855	5-50µL	H,I,J,K,L,M,N,O
W810603	855	20-200µL	I,J,K,L,M,N,O
W810604	855	40-350µL	N,O,P,Q
W810606	855	0.5-10µL	A,B,E,F,G
W810608	855	5-50µL	H,I,J,K,L,M,N,O
W810609	855	20-200µL	I,J,K,L,M,N,O
W810610	855	40-350µL	N,O,P,Q
W851180-2	852	1-10µL	A,B,E,F,G
W870012	926	0.1-2µL	A,B,E,F,G
W870014	926	0.5-10µL	A,B,E,F,G
W870016	926	0.5-10YµL	I,J,K,L,M,N,O
W870018	926	1-20µL	I,J,K,L,M,N,O
W870020	926	2.5-50µL	H,I,J,K,L,M,N,O
W870022	926	5-100µL	H,I,J,K,L,M,N,O
W870024	926	10-200µL	I,J,K,L,M,N,O
W870026	926	50-1000µL	R,S,T,U
W870112	926	0.1-2µL	A,B,E,F,G
W870114	926	0.5-10µL	A,B,E,F,G
W870116	926	0.5-10YµL	I,J,K,L,M,N,O
W870118	926	1-20µL	I,J,K,L,M,N,O
W870120	926	2.5-50µL	H,I,J,K,L,M,N,O
W870122	926	5-100µL	H,I,J,K,L,M,N,O
W870124	926	10-200µL	I,J,K,L,M,N,O
W870126	926	50-1000µL	R,S,T,U
W870128	936	0.1-2mL	V,X,Y
W870130	936	0.25-5mL	X,Y
W870132	936	0.5-10mL	W

Cat. No.	Model	Volume	Tip Ref.
W870134	956	0.5-10µL	A,B,E,F,G
W870136	956	2.5-50µL	H,I,J,K,L,M,N,O
W870138	956	10-200µL	I,J,K,L,M,N,O
W870140	956	20-350µL	N,O,P,Q
W870142	956	0.5-10µL	A,B,E,F,G
W870144	956	2.5-50µL	H,I,J,K,L,M,N,O
W870146	956	10-200µL	I,J,K,L,M,N,O
W870148	956	20-350µL	N,O,P,Q
W870902	926	0.1-2µL	A,B,E,F,G
W870902-A	926	0.1-2µL	A,B,E,F,G
W870904	926	0.5-10µL	A,B,E,F,G
W870904-A	926	0.5-10µL	A,B,E,F,G
W870906	926	0.5-10YµL	I,J,K,L,M,N,O
W870906-A	926	0.5-10YµL	I,J,K,L,M,N,O
W870908	926	1-20µL	I,J,K,L,M,N,O
W870908-A	926	1-20µL	I,J,K,L,M,N,O
W870910	926	2.5-50µL	H,I,J,K,L,M,N,O
W870910-A	926	2.5-50µL	H,I,J,K,L,M,N,O
W870912	926	5-100µL	H,I,J,K,L,M,N,O
W870912-A	926	5-100µL	H,I,J,K,L,M,N,O
W870914	926	10-200µL	I,J,K,L,M,N,O
W870914-A	926	10-200µL	I,J,K,L,M,N,O
W870916	926	50-1000µL	R,S,T,U
W870916-A	926	50-1000µL	R,S,T,U
W870918	936	0.1-2mL	V,X,Y
W870918-A	936	0.1-2mL	V,X,Y
W870920	936	0.25-5mL	X,Y
W870920-A	936	0.25-5mL	X,Y
W870922	936	0.5-10mL	W
W870922-A	936	0.5-10mL	W
W870924	956	0.5-10µL	A,B,E,F,G
W870924-A	956	0.5-10µL	A,B,E,F,G
W870926	956	2.5-50µL	H,I,J,K,L,M,N,O
W870926-A	956	2.5-50µL	H,I,J,K,L,M,N,O
W870928	956	10-200µL	I,J,K,L,M,N,O
W870928-A	956	10-200µL	I,J,K,L,M,N,O
W870930	956	20-350µL	N,O,P,Q
W870930-A	956	20-350µL	N,O,P,Q
W870932	956	0.5-10µL	A,B,E,F,G
W870932-A	956	0.5-10µL	A,B,E,F,G
W870934	956	2.5-50µL	H,I,J,K,L,M,N,O
W870934-A	956	2.5-50µL	H,I,J,K,L,M,N,O
W870936	956	10-200µL	I,J,K,L,M,N,O
W870936-A	956	10-200µL	I,J,K,L,M,N,O
W870938	956	20-350µL	N,O,P,Q
W870938-A	956	20-350µL	N,O,P,Q

¹ 5mL tip used for 2mL unit up to serial number 17091050 ² Does not work on Model 831 ³ Use only on Model 831, will not eject on Model 835 *Graduated