

## EZChrom Elite Compact Hardware and Software Minimum Requirements

The table below lists the hardware and software requirements for EZChrom Elite Compact installations.

- English version of EZChrom Elite Compact is validated on Chinese, Japanese, English, and Western European Language versions of Windows.
- Japanese version of EZChrom Elite Compact is validated on Japanese version of Windows.
- Chinese version of EZChrom Elite Compact is validated on Chinese version of Windows.

| Product                                | Operating System  | CPU (Minimum) | RAM (Minimum) | Network                                    | Disk Space (Minimum)  | Available PC Slots                       |
|--|---|---------------|---------------|--|-----------------------|--|
| EZChrom Elite Compact on Windows XP    | Windows XP Professional (SP3)<br>Microsoft .NET version 3.0 or 3.5<br>SP1 | 2.0 GHz       | 2.0 GB        | Microsoft Networking environment<br>TCP/IP | 10 GB free disk space | Additional Ethernet card may be required |
| EZChrom Elite Compact on Windows Vista | Windows Vista Business (SP1),   | 2.0 GHz       | 2.0 GB        | Microsoft Networking environment<br>TCP/IP | 30GB free disk space  | Additional Ethernet card may be required |

### Important Notes

- EZChrom Elite is validated to run only in a Microsoft Networking environment. Other NOS, such as Unix, Linux or Novell are not supported. WAN's (wide area networks) are not supported.
- Distributed File Sharing (DFS) is not supported for the Enterprise Path or the Enterprise machine.
- RAM is defined as total available memory for EZChrom Elite Compact at any given time. Agilent recommends disabling Windows Vista Aero user interface for Standalone instrument control systems.