


La-Pha-Pack[®]

**HPLC certified
syringes**



La-Pha-Pack HPLC Certified Polypropylene Disposable Syringes

HPLC certified, non sterile disposable syringes made of solvent robust polypropylene for all syringe filter applications

- Each manufactured batch is HPLC controlled and supplied with an appropriate certificate upon request. The certificate clearly states column and other method conditions being used for the test.
- Syringes manufactured according to ISO 13485 (Medical devices)
- Luer Lock manufactured according to ISO 594-2 / DIN EN 1707
- Luer Slip manufactured according to ISO 594-1 / DIN EN 20594-1
- Two-part, all-plastic construction made of a chemically resistant, inert polypropylene
- No rubber plunger seals or silicone lubricants that may cause sample contamination
- Free of latex, free of plasticizers, free of PVC
- Safe functional back-stop feature
- Easy to read permanent graduations
- All Luer-Lock syringes have centered tips
- Available with Luer-Lock and Luer slip connections

See back page for ordering information 

Combine the syringes with our high-quality HPLC syringe filters.

Available in sizes 17mm, 25mm and 30mm. Further cost-free test samples upon request!

17mm



25mm



30mm



Order information:

HPLC certified
syringes



02 36 2399



05 36 2400



10 36 2401



20 36 2402

Luer-Lock

La-Pha-Pack Disposable Syringes with Luer-Lock

Description	Volume (ml)	Part Number	Packaging
Disposable syringes, non sterile, Luer-Lock	2ml	02 36 2399	100 pcs. per PE-Bag
Disposable syringes, non sterile, Luer-Lock	5ml	05 36 2400	100 pcs. per PE-Bag
Disposable syringes, non sterile, Luer-Lock	10ml	10 36 2401	100 pcs. per PE-Bag
Disposable syringes, non sterile, Luer-Lock	20ml	20 36 2402	100 pcs. per PE-Bag

HPLC certified
syringes



02 36 2403



05 36 2404



10 36 2405



20 36 2406

Luer-Slip

La-Pha-Pack Disposable Syringes with Luer-Slip

Description	Volume (ml)	Part Number	Packaging
Disposable syringes, non sterile, Luer-Slip	2ml	02 36 2403	100 pcs. per PE-Bag
Disposable syringes, non sterile, Luer-Slip	5ml	05 36 2404	100 pcs. per PE-Bag
Disposable syringes, non sterile, Luer-Slip	10ml	10 36 2405	100 pcs. per PE-Bag
Disposable syringes, non sterile, Luer-Slip	20ml	20 36 2406	100 pcs. per PE-Bag

For more information, please contact customer service:

Tel.: +49 2423 94 31 – 0

Fax: +49 2423 94 31 – 34

e-mail: service@la-pha-pack.com

www.la-pha-pack.com