

KIT CONTENT

| Order No. | Volume | Symbol | Component |
|---------------------|------------|----------------|---|
| ZV-3046-02S1-10 | 1 x 2 ml | STD I | Serum Calibrator, lyophilised |
| ZV-3046-02K1-10 | 1 x 2 ml | CONT I | Serum Control, lyophilised |
| ZV-3046-02K2-10 | 1 x 2 ml | CONT II | Serum Control, lyophilised |
| ZV-3046-02R1-10 | 2 x 50 ml | R1 | Reagent 1, Contains aqueous buffer |
| ZV-3046-02R2-10 | 2 x 50 ml | R2 | Reagent 2, Contains deuterated internal standard |
| ZV-3046-02MP-10 | 1 x 500 ml | MP | Mobile Phase, Contains organic solvent |
| ZV-3046-02MT-10 | 1 x 500 ml | SPE | SPE buffer Contains organic solvent |
| ZV-3046-02WB-10 | 1 x 500 ml | WB | Washing Solution Contains organic solvent |
| ZV-3046-KK-En_Rev01 | 1 x 1 pc | KK | User Guide |