

## KIT CONTENT

Order No.	Volume	Symbol	Component
ZV-3017-02K1-15	1 x 1 Spot	<b>CONT I</b>	Dried Blood Control I, Contains Human Blood
ZV-3017-02K2-15	1 x 1 Spot	<b>CONT II</b>	Dried Blood Control II, Contains Human Blood
ZV-3017-05R1-15	2 x 50 ml	<b>R1</b>	Reagent 1, Contains deuterated internal standards
ZV-3017-05R2-15	1 x 30 ml	<b>R2</b>	Reagent 2, Contains derivatisation reagent
ZV-3017-05R3-15	2 x 50 ml	<b>R3</b>	Reagent 3, Contains acid
ZV-3017-05WB-15	1 x 500 ml	<b>WB</b>	Washing solution, Contains organic solvent
ZV-3017-05MP-15	1 x 1000 ml	<b>MP</b>	Mobile Phase, Contains organic solvent
ZV-3002-05KK15	1 x 1 pcs.	<b>KK</b>	User Guide