

KIT CONTENT

Order No.	Volume	Symbol	Component
ZV-3046-02S1-10	1 x 2 ml	STD I	Serum Calibrator, lyophilised
ZV-3046-02K1-10	1 x 2 ml	CONT I	Serum Control, lyophilised
ZV-3046-02K2-10	1 x 2 ml	CONT II	Serum Control, lyophilised
ZV-3046-02R1-10	2 x 50 ml	R1	Reagent 1, Contains aqueous buffer
ZV-3046-02R2-10	2 x 50 ml	R2	Reagent 2, Contains deuterated internal standard
ZV-3046-02MP-10	1 x 500 ml	MP	Mobile Phase, Contains organic solvent
ZV-3046-02MT-10	1 x 500 ml	SPE	SPE buffer Contains organic solvent
ZV-3046-02WB-10	1 x 500 ml	WB	Washing Solution Contains organic solvent
ZV-3046-KK-En_Rev01	1 x 1 pc	KK	User Guide