

Item No.	Description
----------	-------------

Glassware

14-2052-024	Sparger, 5mL needle, 1/2" (glassware only)
14-2053-024	Sparger, 25mL needle, 1/2" (glassware only)
14-2333-024	Sparger, 25mL fritless, 1/2" (glassware only)
14-2334-024	Frit Sparger, 25mL, 1/2" neck, Tekmar
14-2336-024	Fritless Sparger, 5mL, 1/2" Tekmar
14-2337-024	Frit Sparger, 5mL, 1/2" Neck Tekmar
14-8921-000	Sparger kit, 5mL fritted
14-8922-000	Sparger kit, 25 mL fritted
14-8923-000	Sparger kit, 5 mL fritless
14-8924-000	Sparger kit, 25 mL fritless
14-8925-100	kit, 5mL needle sparge, Stratum
14-8926-100	kit, 25mL needle sparge, Stratum

Heaters

14-7806-000	Assy, RTD (24") Temp Sensor for Soil Valve, Mount
14-8785-000	Oven RTD assembly
14-8788-300	Assy, Oven Heater, 115V
14-8789-200	Assy, Mount HTR 115V
14-8789-300	Mount Cartridge Heater 230V (Stratum)
14-8790-C20	Transfer Line Heater Assy, 72", 115V
14-8790-D20	Transfer Line Heater Assy, 72", 230V
14-8920-200	Sample Heater Assembly for all glassware types, 110V
14-8788-300	Assy, Oven Heater, 115V
14-9828-220	Trap Heater, Nichrome Wire, 230V
14-9895-120	Heater, Condensate Trap, 115V
14-9895-220	Heater, Condensate Trap, 230V
14-8920-300	Assy, Heater, Sample 5/25mL 230VAC
14-9828-120	Trap Heater, Nichrome Wire, 115V

Fittings

12-0041-016	Ferrule, teflon, 1/4", glassware
12-0044-016	Ferrule set, 1/8", brass, set
12-0069-016	Nut, 1/8", brass
12-0070-016	Tee, 1/8", Brass
14-0051-016	Union, 1/16-1/16 SS
14-0053-016	Elbow, male, 1/8", NPT, brass
14-0158-016	Ferrule, 1/16", SS, Set
14-0216-016	Female luer connector for sample valve
14-0241-016	Ferrule, 1/16", SS
14-0243-016	Nut, 1/16", short, SS, Valco
14-1590-116	Nut, Male Plug, 1/16" Gold Plated
14-2261-016	Union, 1/4 -1/16 Reducing, SS, W/ Teflon Fer
14-2549-016	Elbow, 90 degree, tube, 1/8 NPT to 1/16
14-3295-016	Nut, 1/16", male, SS, Swagelok
14-4050-016	Ferrule, 1/8 Fitt 1/8 Tube GVF
14-8470-016	Plug Nut, 1/4-28, Tefzel, Flat Bottom
14-8914-016	Nut, 1/8 Black Delrin Upchurch
14-8915-016	1/8" Manifold Bushing, Delrin
14-3098-016	Ferrule, 1/2", Teflon, Valco
14-9801-016	Mount/Tee Assy, Siltek w/ Fittings (no heaters)
14-5773-016	Ferrule, 1/8", Gold Plated for Top of Trap
14-9845-016	Trap Manifold Assy w/ Fittings
885-291	Check valve 1/8"

14-4602-116	Nut, 1/8", SS, Valco Gold Plated
14-8820-016	Nut, 1/16", Short, PEEK
14-8659-016	Nut, Flangless, Short, 1/8, Peek
14-0243-116	Nut, 1/16", Short Gold plated
14-7671-016	Ferrule, Super Flangeless, 1/16", PEEK
14-7695-016	Nut, PEEK, Natural for 1/16" OD Tubing
15-0285-016	1/8" Vespel ferrule
14-5661-016	Pencil Filter, Filter Flow Module

Electronics

080-613	Battery, Lithium, SIZE 2032
14-7835-190	Temperature Control Board
14-8781-190	PCB, DC Control
14-5665-034	Fuse, Slow Blow, 5AMP, 250V, 5 X 20 MM
14-7837-590	Microprocessor (CPU) Board for Stratum
14-9185-090	PCBA, AC Output Board
14-8312-090	Power Supply, triple out, +5/+15/-15 VDC
14-9577-090	PCB, LED Banner Control Assembly
14-4383-028	Switch, Power, 10 AMP Filtered use Delta
14-7682-051	Solid State Relay, 240 VAC, 25A, 0 Volt
14-8834-155	Foam Module, All Units, Green
14-9840-100	Trap Blower Fan Assy
14-9852-090	Power Supply, +24V, 65W

Valves

14-0036-050	Sample valve, 3-port
14-6154-350	6 port control box only
14-6292-350	Controller Box, Valve, 6 Port, Velocity 6-port Valve Assembly, small body, 24 VDC includes Valve, fittings, actuator, controller and cable
14-9126-050	6-port Valve small body w/ slider.
14-9126-150	6 port actuator only, SM body valve
14-9126-250	Rotor, 6-port valve, small body after SN US04119002
14-9126-350	Drain/Foam 24VDC Solenoid Valve Kit
15-0485-000	3-port solenoid valve for manifold
14-8901-750	Mass Flow Controller Assy
14-8906-050	Solenoid Valve Manifold Assy, LQX

Miscellaneous

14-1197-039	Cord, Power, 115V 15 A, European Recept
14-2440-006	Ceramic Standoff for valve oven plate
14-8856-000	Assy, 10mL Antifoam 1520
14-8909-043	Thermal Green Septa for Trap Injection Port
14-9847-000	Cooling Fan Assy, Low CFM
14-8834-000	Guardian Foam Sensor option
14-8834-100	Guardian Foam Sensor and Eliminator Retrofit Kit VOC Teklink Control Software capable of controlling the StratUm PTC with
14-9800-075	AQUATek 70 / SOLATek 72 autosampler systems
14-9800-PRT	StratUm PTC Refurb Kit, parts only
14-9800-RFB	StratUm PTC Refurb Kit, with labor
14-9809-054	65lb Gas Shock
14-9800-PTLV	Refurb Kit, Stratum, Parts Only, W/O 6-Port
14-9849-000	Top of Trap Injection Port Assy (StratUm)
14-5028-039	220V Power Cable
15-0439-000	External Pressure Regulator
14-9800-KIT	Kit box, Stratum

Traps

12-0083-403	Trap, U Shape, Tenax 12" LG, Stamp "1"
12-0084-403	Trap, Tenax/SG, 12" LG, Stamp "2", U Shape
14-0124-403	Trap, Tenax/SG/Char 1/8" Stamp "3", U Shape
14-2366-403	Trap, Char-SC-TNX-OV1 12X1/8, "#5" U Shape
14-3347-403	Trap, OV-1. Tenax, #7 U Shape
14-3928-403	Strat-Trap, Carbo B/Carbo S, (#8)
14-5865-403	Strat-Trap, Vocarb 4000
14-5866-403	Strat-Trap, BTEX
14-9333-403	Strat-Trap, BTEX + MTBE
14-9908-403	Strat-Trap, #9 (U-shaped)
14-9909-403	Strat-Trap, Tenax/Silica Gel/Carbo Sieve III (#10)
14-5864-403	Strat-Trap, Vocarb 3000

Tubing

14-6934-002	PEEK Tubing, 1/16" OD, Natural (per foot)
14-7630-002	Tubing, 1/8" OD, Red, Translucent (per inch)
14-7634-002	Tubing, 1/8" OD, Blue, Translucent (per inch)
14-8927-002	PEEK Tubing, 1/16" OD, Green (per foot)
14-8928-002	PEEK Tubing, 1/16" OD, Orange (per foot)
14-9622-002	Tubing, Siltek, 1/16" OD x 0.040" ID (per foot)
14-9623-002	Tubing, Transfer Line, Siltek (per foot), 0.021 ID x .029 OD
14-9927-002	Tubing, Siltek, 1/16" OD x 0.02" ID (per foot)

GC Interface Cables

14-2991-000	Interface, HP 5890 GC
14-2993-000	Interface, HP 5995/96/85/87/88/92 GC/MS with HP-1000/RTE GC/MS Software.
14-3147-000	Interface, general purpose/HNU 301/321/421 (Valve driver option necessary from HNU.)
14-3233-000	Interface, Perkin-Elmer 8000 Series/Autosystem, Clarus.
	Interface, HP 5890 w/5970 MSD and Unix or Pascal based software. Note: This cable requires the HP A111 (HP P/N 05990-60111) or A211 (HP P/N 05990-60211) accessory card and internal accessory cable (HP P/N 05987-60158) if using the Pascal software. If using Unix software the internal accessory cable (HP P/N 05987-60158) only is necessary.
14-4188-086	Interface, Shimadzu GC 14A/15A/17A GC 14A w/QP 1000 EX MSD and GC 14A w/QP 2000 MSD, also includes 1/16" to 2mm unions, needs Shimadzu P/N 221-34280-91
14-4610-086	Interface, Shimadzu GC 14A/15A/17A GC 14A w/QP 1000 EX MSD and GC 14A w/QP 2000 MSD, needs Shimadzu P/N 221-34280-91
14-4610-186	Interface, HP 5890 w/5970 MSD and Unix-B or MS-DOS software, HP 5890/5971/5972 MSD and Unix-B or MS-DOS software and HP 5890/5989 MS Engine
14-4652-086	Interface, Carlo Erba Mega and Vega Series/Fisons 8000
14-4938-086	Interface, two Tekmars on one Perkin-Elmer 8000 Series/Autosystem
14-5397-086	Interface, HP 6850, 6890, 7890 GC
14-6689-086	Interface, HP 6850, 6890, 7890 GC
14-7396-086	Interface, Varian 3800/3900
14-7397-086	Trace 2000/Trace GC Ultra Interface cable w/ 1/16" to 2mm unions
14-7397-186	Trace 2000/Trace GC Ultra Interface cable
	Duet™ High Throughput Interface. Interfaces two Teledyne Tekmar VOC systems with any GC or GC/MS system.
14-8417-000	PE Clarus 500 Interface Cable
14-3233-200	PE Clarus 500 Interface Cable
14-9754-086	Shimadzu 2010-2014 Interface Cable
14-8477-286	Cable, Archon-Tekmar Interface

Cables



14-8876-086	Mass Flow Controller Cable
14-8797-286	Cable Assy, CPU SPI Interface, StratUm
14-8798-286	Cable Assy, AA70 to Concentrator
14-8796-186	Short LNET Cable
14-9857-086	DC power harness cable assy
14-9858-086	AC power harness cable assy
14-8020-186	Cable Assy, Temp to SS relay ribbon
14-9880-086	Trap Heater RTD extension cable assy
14-8901-550	Valve LED Power Cable, W/3" LED Ext.
14-8901-450	LED Cable, Drain Valve
14-9654-086	LED Banner cable assy
14-9879-086	Trap AC extension cable, 115V