



# Replacement Parts for the Agilent Dual Plasma Burner and Controller for Use with the 355 SCD and 255 NCD

Part Number	Part Description	Unit Code
Burner Components		
G6600-60037	Dual Plasma Burner Accessory Kit, SCD	EA
G6600-60038	Dual Plasma Burner Accessory Kit, NCD	EA
G6600-80018	Dual Plasma Burner Column Nut and Ferrule	EA
G6600-80046	DP Valcon Ferrule 1/32 x .5mm	EA
G6600-80066	DP Graphite 1/4 inch Ferrule	EA
G6600-60614	DP 6890 Heater/Sensor (for Agilent 6890)	EA
G6600-60014	Dual Plasma Thermocouple Probe, K-Type	EA
G6600-20000	Dual Plasma Burner Splitter Fitting	EA
G6600-20007	Fitting DP Burner Inlet Assembly	EA
G6600-20004	DP Tapered Union Fitting	EA
G6600-60006	Quartz Burner Heater Assembly	EA
G6602-64030	Complete Spare SCD Dual Plasma Burner w/tubes	EA
G6600-64030	Complete Spare NCD Dual Plasma Burner w/tubes	EA
Controller Components		
G6600-20000	SS Gas Delivery tube for DP	EA
G6600-60013	DP Controller Internal Tubing Kit	EA
G6600-60036	Dual Plasma Controller Solenoid Valve	EA
G6600-80049	DP Controller Relay, Solid State	EA
G6600-67000	DP Controller CPU Replacement Kit	EA
G6600-60056	DP Controller Flow Controller	EA
G6600-67001	DP Transformer Replacement Kit	EA
G6600-60012	DP Controller In/Out Hardware Kit	EA



<b>Part Number</b>	<b>Part Description</b>	<b>Unit Code</b>
Detector Components		
G6600-67004	SCD UV Pass Filter Replacement Kit	EA
G6600-67003	NCD Red Pass Filter Replacement Kit	EA
G6600-60051	Unheated Transfer Line, 1/8 Inch	EA
G6600-80033	Ozone Post Restrictor	EA
G6600-67002	Ozone Generator Replacement Kit	EA
G6600-60000	Ozone HV Transformer – 100/120V HV	EA
G6600-60001	Ozone HV Transformer - 230V HV	EA
Pumps & Components		
G6600-85000	Chemical Trap	EA
G6600-80044	Replacement Oil Coalescing Filter	EA
G6600-85001	Synthetic Oil for Edwards Pump	EA
G6600-60021	Repair Kit for Dry Piston Pump	EA
G6600-64020	Vacuum Pump,100VAC, Dry Piston	EA
G6600-64021	Vacuum Pump,115VAC, Dry Piston	EA
G6600-64022	Vacuum Pump, 230VAC, Dry Piston	EA
Other		
G2933-85003	Sulfur Trap	EA
G6600-60600	DP Burner 6890 FID Mounting Kit Adapter	EA
G6600-60623	DP Burner 5890 FID Mounting Kit Adapter	EA
G6600-67007	Maintenance Kit, 12 Month Dual Piston RV5 Pump	EA
G6600-67008	Maintenance Kit, 12 Month Dry Piston Pump	EA



Copyright © 2007  
Agilent Technologies, Inc.  
2850 Centerville Road  
Wilmington, DE 19808-1610  
Printed in the USA Feb 2007



G6600-90030