

La-Pha-Pack®



La-Pha-Pack Total Phthalate Free Seal (T P F Seal) with assembled Aluminum Liner

La-Pha-Pack offer the first 11mm Crimp Cap with assembled aluminum liner

- Septa material free of any elastomers and halogens; therefore no contamination of the sample by plasticizers, Silicone or Butyl rubbers and PTFE / FEP / TEF materials
- Recommended for GC/MS and LC/MS applications
- The additional ring which is assembled above the aluminum liner achieves an excellent tightness and inertness of the complete combination seal
- Suitable for storing of standards or reactive substances which would normally corrode the septum material
- Packed with 100 pieces in a packaging free of plasticizers (glass container)

Application areas:

- Analytic of elastomers und plastomers
- Phthalate analysis
- Analytic of fluorinated, halogenated organic compounds
- VOC (volatile organic compounds) analysis
- Analytic of polymerisation catalysts
- Analytic of silicones und silicates



Ordering information:

Article No.: 11 03 2519

11mm TPF Combination Seal: Aluminum Cap, clear lacquered, centre hole; aluminum liner, 0.06mm (Sealing by an additional installed transparent silicone-flat seal)

For more information, please contact our customer service:

Tel.: +49 2423 94 31 – 0

Fax: +49 2423 94 31 – 34

e-mail: service@la-pha-pack.com

www.la-pha-pack.com

©2012 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

Part of Thermo Fisher Scientific