

UNIS 3100 HT PTV



The perfect choice for high temperature applications

- S/SL and PTV modes
- max. temperature: 500 °C
- heating rate: 720 °C/min
- low thermal mass allows rapid heating and cooling
- extended heater to avoid 'cold spots' - even at an oven temperature of -20 °C

This UNIS inlet is designed for special high temperature applications like SimDis or H53. The special design with its low thermal mass is based on the UNIS 2100 PTV injector and can reach a maximum temperature of 500 °C. The UNIS 3100 HT PTV provides both split/splitless and PTV modes.