



La-Pha-Pack HPLC and GC Certified Vial KITS

Certifications are getting more and more important, in order to make processes more reproducible and avoid possible sources of errors from the beginning. For La-Pha-Pack highest quality, consistency and quality control have always been very important and have been documented by HPLC and GC certified vials and closures.

- Each batch of HPLC and GC-certified KITS is tested on 15 critical parameters. In a realistic method, a HPLC/UV and GC/FID test of vials/closure combination for blanks and contaminations will be carried out.
- The batch-specific test certificate with the HPLC and GC-Chromatograms can be handed out upon request.
- The HPLC and GC certified KITS are delivered completely shrinkwrapped for reasons of originality, purity and transport safety. This means an additional safety for the end user.
- Available as 9mm Short Thread Vial in clear and amber with suitable closure combination.
- Upon request further HPLC and GC Certified Vial KITS will be available.

Ordering information

La-Pha-Pack Certified Vials Kits

Description	Article No.
HPLC/GC Certified 1.5ml Short Thread Vial, clear glass, 1st hydrol. class, with label; UltraClean Seal: 9mm PP Short Thread Cap, blue, centre hole; Silicone white/PTFE red, 1.0mm	11 40 2556
HPLC/GC Certified Vial KIT 1.5ml Short Thread Vial, amber glass, 1 st hydrol. class, with label; UltraClean Seal: 9mm PP Short Thread Cap, blue, centre hole; Silicone white/PTFE red, 1.0mm	11 40 2557

For more information, please contact customer service:

Tel.: +49 2423 94 31 – 0

Fax: +49 2423 94 31 – 34

e-mail: service@la-pha-pack.com

www.la-pha-pack.com

©2012 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.