


Agilent 1260 Infinity Manual Injector

Features and Ordering Details



For low bench space footprint

Convenient modular accessory, dimensioned to match the footprint of all Agilent 1260 Infinity modules. Stack to create a compact, affordable manual injection liquid chromatography. It is compatible with up to 600-bar operating pressure.

Features

- Lowest bench space requirements and highest flexibility with external HPLC columns by using the integrated column and base
- Long valve life, up to 600-bar operating pressure
- High instrument uptime through proven parts, e.g. ceramic rotor faces



Agilent Technologies

Ordering Details – Agilent 1260 Infinity Manual Injector

Description	Product Number
Agilent 1260 Infinity manual injector with 20 μ L sample loop and mounting stand	G1328C
5- μ L sample loop	0101-1248
10- μ L sample loop	0100-1923
20- μ L sample loop	0100-1922
50- μ L sample loop	0100-1924
100- μ L sample loop	0100-1921
200- μ L sample loop	0101-1247
500- μ L sample loop	0101-1246
1-mL sample loop	0101-1245
2-mL sample loop	0101-1244
5-mL sample loop	0101-1243

www.agilent.com/chem/1200

© Agilent Technologies, Inc., 2010-2012
Published in USA, November 1, 2012
Publication Number 5990-6122EN



Agilent Technologies