

La-Pha-Pack®



La-Pha-Pack PP Storage Boxes


La-Pha-Pack PP Storage Boxes for a convenient and safe storage and transportation of sample vials

- Safe standing position on the laboratory table and during transport due to specific cavities related to the vial diameter
- Ideal for space-saving storage in fridges, as the transparent lid prevents condensations on the closures and thus avoids a possible contamination in the cooling unit.
- Temperature resistant from -80 degrees / + 100 degrees
- Alpha Numeric Coding (1.5 ml, 4ml) for clear sample identification
- Unbreakable polypropylene bottom and lid, stackable
- Chemically resistant and fairly robust; autoclavable

Available for:



**1.5ml (1.8ml, 2ml) Crimp, Snap and Screw Neck Vials (ND8, ND9, ND10, ND11)
2ml Shell Vials**



4ml Crimp and Screw Neck Vials (ND13)



5ml, 10 ml, 20ml Headspace Vials (ND18 / ND20)



20ml, 30ml, 40ml EPA Vials (ND24)

See back page for ordering information 

For more information, please contact our customer service:

Tel.: +49 2423 94 31 – 0

Fax: +49 2423 94 31 – 34

e-mail: service@la-pha-pack.com

www.la-pha-pack.com

©2012 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

Part of Thermo Fisher Scientific

Ordering information:




PP Storage Boxes for 1.5ml Sample Vials

Description

PP Storage Box for 1.5ml (1.8ml,2ml) vials or 2ml shell vials, blue, with cover, (130 x 130 x 45mm) , 81 cavities with alpha numeric coding of all 4 margins as well as the cavities at the bottom

Article No.

BLUE	ORANGE	PINK	YELLOW	GREEN	TRANSP.
12 21 2420	12 21 2421	12 21 2422	12 21 2423	12 21 2424	12 21 2425



PP Storage Box for 4ml Sample Vials

Description

PP Storage Box for 4ml vials or 4ml shell vials, red, with cover, (130 x 130 x 52mm) , 49 cavities with alphanumeric coding at the cavities

Article No.

15 21 2426



PP Storage Boxes for 5ml, 10ml and 20ml Headspace Sample Vials

Description

PP Storage Box for 5ml, 10ml and 20ml Headspace - vials, blue, with cover, (130 x 130 x 102mm), 25 cavities

Article No.

24 21 2427



PP Storage Box for 20ml, 30ml, 40ml EPA - Vials, violet, with cover

Description

PP Storage Box for 20ml EPA - vials, violet, with cover, (130 x 130 x 102mm), 16 cavities

PP Storage Box for 30ml and 40ml EPA-Vials, violet, with cover (130 x 130 x 105mm), 10 cavities

Article No.

28 21 2428

28 21 2429