



## La-Pha-Pack Pink Caps

### La-Pha-Pack offer new "PINK" Autosampler Vial Caps

- Available as 8mm and 9mm screw thread and 11mm snap caps.
- Reliable and trusted options for any autosampler and any application.
- Vial and cap are also available as 2in1 KITS.
- Caps are available with various liner materials, such as Natural Rubber/Teflon, RedRubber/PTFE, Silicone/PTFE, Silicone/PTFE slitted and PTFE/Silicone/PTFE.

#### Available as:



#### 8mm Screw Neck Seals ND 8, open top

Article No. 08 15 2466	Article No. 08 15 2461	Article No. 08 15 2463	Article No. 08 15 2464
Natural Rubber red-orange/TEF transp., 60° shore A, 1.3mm	RedRubber / PTFE beige, 45° shore A, 1.0mm	Silicone white/PTFE red, 45° shore A, 1.3mm	Silicone cream/PTFE red, 55° shore A, 1.5mm



#### 9mm Short Thread Cap ND 9, open top

Article No. 09 15 2469	Article No. 09 15 2470	Article No. 09 15 2471	Article No. 09 15 2472
Natural Rubber red-orange/TEF transp., 60° shore A, 1.3mm	RedRubber / PTFE beige, 45° shore A, 1.0mm	PTFE red/Silicone white/PTFE red, 45° shore A, 1.0mm	Silicone white/PTFE blue, 55° shore A, 1.0mm, slitted



#### 11mm Snap Ring Cap ND 11, open top

Article No. 11 15 2474	Article No. 11 15 2478	Article No. 11 15 2475	Article No. 11 15 2477
Natural Rubber red-orange/TEF transparent, 60° shore A, 1.0mm	RedRubber / PTFE beige, 45° shore A, 1.0mm	PTFE red/Silicone white/PTFE red, 45° shore A, 1.0mm	Silicone white/PTFE blue, 55° shore A, 1.0mm, slitted

For more information, please contact our customer service:

Tel.: +49 2423 94 31 – 0

Fax: +49 2423 94 31 – 34

e-mail: [service@la-pha-pack.com](mailto:service@la-pha-pack.com) / [www.la-pha-pack.com](http://www.la-pha-pack.com)