

**New:**



## FID for Agilent 1260 Infinity SFC

- Fast analysis of fuels according to ASTM D5186
- Seamless integration in OpenLAB CDS ChemStation
- Ease of use
- Low operating costs by using standard grade gaseous CO<sub>2</sub>
- Green chemistry: low solvent consumption and less waste generation

Distribuido en España por:


 **Ingeniería Analítica**  
Chromatography & Spectrometry

Tel: (+34) 93.590.28.50 Fax: (+34) 93.675.05.16 [www.ingenieria-analitica.com](http://www.ingenieria-analitica.com)  
[inf@ingenieria-analitica.com](mailto:inf@ingenieria-analitica.com)


**SIM**   
Scientific Instruments Manufacturer GmbH

[www.sim-gmbh.de](http://www.sim-gmbh.de) phone +49-208-9410780

 **Agilent Technologies**  
Premier Solution Partner



D5186 test mix



Screenshot OpenLAB CDS ChemStation